



Material Safety Data Sheet: Tuffie Detergent Wipes

1.0 Identification of Preparation AND Supplier

Product Name: Tuffie Detergent Wipes (225 Tub)
Product Code: 901DW225
Recommended Use: Multi-surface wipes saturated with a low foaming, amphoteric cleaning solution for cleaning a number of dirty and dusty surfaces.
Review / Issue No / Date 3 / October 2009

Manufacturer: **Vernacare Ltd**
Address: Folds Rd, Bolton, BL1 2TX
Telephone No: + 01204 529494

Australian Importer: **EBOS Healthcare**
Address: Unit 2, 109 Vanessa St, KINGSGROVE, NSW 2209
Telephone No: Free Call 1800 269 534

2.0 Hazards Identification

Classification: This product is not regarded as a health hazard under current legislation.
Health Hazards: Contains components classified as Irritant
Environmental Hazards: Not classified as dangerous under EC criteria.
The detergent components of this product are biodegradable
Safety: This material is non combustible.



Material Safety Data Sheet: Tuffie Detergent Wipes

3.0 Composition/Information on ingredients

Wipes: Thermo bonded blend of polypropylene and viscose fibres.

Preparation: Irradiated deionised water, Ethylan CPG, Betain 810, Liquid Glydant Plus, 38% Tetra sodium EDTA solution, Citric Acid. Contains less than 5% non-ionic and amphoteric surfactants preservatives.

4.0 First Aid Measures

Eye Contact: Avoid eye contact. If contact occurs, Irrigate thoroughly with water for at least 15 minutes. If persistent irritation occurs, obtain medical attention.

Skin Contact: Avoid prolonged contact. Wash hands thoroughly with water after use.

Ingestion: May cause irritation to mucous membranes. Do not induce vomiting. If irritation or symptoms persist, seek medical attention

5.0 Fire Fighting Measures

Extinguishing Media: Use any available extinguisher.

Special Exposure Hazards: None known

Protective Equipment: Wear full fire fighting protective clothing including self contained breathing apparatus

6.0 Accidental Release Measures

Personal Precautions: Avoid contact with eyes

Environmental Precautions: Dispose of wipes in suitable containers

Clean Up: This product cannot be spilled



Material Safety Data Sheet: Tuffie Detergent Wipes

7.0 Handling and Storage

- Handling:** No special handling requirements.
- Storage:** Protect from extremes of temperature and store out of direct sunlight
- Shelf Life:** Three years from date of delivery

8.0 Exposure Controls/Personal Protection

- Exposure limit values:** None
- Exposure controls:** None
- Respiratory Protection:** Not required in normal use.
- Hand Protection:** Not required under normal conditions of use.
May cause irritation to sensitive skin.
Gloves advisable for sensitive skin.
- Eye Protection:** Not required in normal use.
- Skin Protection:** Not required under normal conditions of use.
May cause irritation to sensitive skin.

9.0 Physical & Chemical Properties

- Appearance:** White, moist and soft containing a liquid free from matter
- Wipe Dimensions:** 300mm x 250mm
- Wipe Weight:** 50gsm
- Odour:** Odourless
- pH:** 5.5
- Miscibility:** Miscible in water



Material Safety Data Sheet: Tuffie Detergent Wipes

Flash Point: N/A

10.0 Stability and Reactivity

Stability: Stable

Conditions to Avoid: Extremes of temperature. Direct sunlight.

Materials to Avoid: Do not mix with any other cleaning chemicals

Decomposition Products: N/A

11.0 Toxicological Information

Classification: Not classified – no specific health warning noted.

Eye Contact: May cause irritation.

Skin Contact: May cause skin dryness and cracking with prolonged contact

Inhalation: Not a likely route of exposure.

Ingestion: May cause irritation to the mouth and upper digestive tract.

12.0 Ecological Information

Biodegradability: The product is inherently biodegradable. No environmental hazard is anticipated with small volumes associated with this product, provided that the product is handled and disposed of with due care and attention.

13.0 Disposal Considerations

No specific disposal methods are required but the wipes should not be flushed away.

Dispose of waste material in waste container



Material Safety Data Sheet: Tuffie Detergent Wipes

14.0 Transport Information

Hazard Classification: Not classified, no warnings

15.0 Regulatory Information

Classification: Not classified

EC Classification: Not classified

EC Risk Phrases: Not classified

UK Regulatory References: Health and Safety at Work Act 1974

Additional Safety Information: As with all detergents, people with sensitive or damaged skin should avoid prolonged contact with the product

Country of Origin: United Kingdom

16.0 Other Information

Training: It is recommended that all personnel are trained in the safe handling of chemicals.

Recommended Use: Wipes for cleaning Surfaces and Equipment

Product Description: 901DW225 Tuffie Cleaning Wipe

The information contained in this safety data sheet is given in good faith. It is accurate to the best of our knowledge and belief and represents the most up to date information.

The information given in this data sheet does not constitute or replace the users own assessment of workplace risk as required by other health and safety legislation.