



Tuffie

cleaning and disinfecting wipes
the complete range



Vernacare

 **EBOS**
HEALTHCARE

ESSENTIAL TO HEALTH




Tuffie detergent wipes

A range of wipes for effective cleaning of hard surfaces and general equipment.

Perfect for using where general cleanliness is required quickly and conveniently.

Specially formulated to remove dirt from surfaces and equipment:



Ideal for use on mattresses, beds, commodes, door handles, hoists, trolleys and other general equipment.

Other benefits when using Tuffie Detergent Wipes:

- Improved cleaning performance
- Non-streak formula
- Faster drying time
- Active Odour Control
 - Attacks & destroys bad smells
 - Effective against vomit, faeces and urine
 - No allergen labelling required

1. Tuffie Detergent Wipes - Bucket of 225

EBOS Code: 10020672

Supplier Code: 901DW225

Bracket EBOS Code: 10009577

Supplier Code: TTWB225/2275

2. Tuffie Detergent Wipes Flexican - 6 x 150's

EBOS Code: 10020671

Supplier Code: 901DW150RP

Bracket/Holder

EBOS Code: 10020670

Supplier Code: 901DW150FC

3. Tuffie Detergent Wipes Refill - 6 x 150's

EBOS Code: 10020671

Supplier Code: 901DW150RP

4. 'Dream' Dispenser

EBOS Code: 10020760

Supplier Code: 520ZS001

Vernacare

EBOS
HEALTHCARE



Images for illustrative purposes only.

**Alcohol
FREE**

**Tuffie 5 is
effective on
Coronavirus**

Tuffie 5

CLEANING & SANITISING WIPES

LOW LEVEL INSTRUMENT GRADE DISINFECTANT

Surface cleaning and disinfection is one of the most effective ways to stop the spread of infection. Various studies have shown that contaminated surfaces are an established route of transmission for high-risk pathogens.¹ The WHO advise that standard precautions for all patients are applied in the areas of waste management and environmental cleaning.

Tuffie Cleaning and Sanitising Wipes:

- Powerful alcohol and chlorine-free wipes, ideal for cleaning and disinfecting all non-porous surfaces
- Tested to EN standards
- Quick drying no residue formulation helps reduce the risk of cross-infection within busy healthcare environments.

**Cleans and sanitises
all hard surfaces.**



**Alcohol
Free**



**Fragrance
Free**



**Standard
Approved**

Micro-Organism	Standard
Acinetobacter baumannii (Multidrug-resistant) ATCC BAA-1605, Enterococcus faecalis ATCC 19433, Enterococcus faecalis (Vancomycin-resistant) ATCC 51299, Enterococcus faecalis (Vancomycin-resistant) ATCC 51575, Enterococcus hirae ATCC 10541, Escherichia coli ATCC 10536 Klebsiella pneumoniae (ESBL) ATCC 700603 Listeria monocytogenes ATCC 13932, Pseudomonas aeruginosa ATCC 15442, Staphylococcus aureus ATCC 6538, Staphylococcus aureus (Methicillin-resistant) ATCC 33592	Bactericidal - EN 13727
Pseudomonas aeruginosa ATCC 15442, Staphylococcus aureus ATCC 6538, Enterococcus hirae ATCC 10541	Bactericidal - EN 16615
Mycobacterium terrae ATCC 15755	Mycobactericidal - EN 14348
Candida albicans ATCC 10231	Bactericidal - EN 16615
Candida albicans ATCC 10231	Fungicidal - EN 13624
Candida auris NCPF 8984	Fungicidal - EN13624
Feline coronavirus, strain NADL-2, FLI-RVB-1259	Virucidal - EN 14476
Vaccinia virus, strain Ankara, ATCC VR-1508	Virucidal - EN 14476
Human coronavirus, strain 229E, ATCC VR-740	Virucidal - EN 14476
Adenovirus type 5, strain Adenoid 75, ATCC VR-5	Virucidal - EN 14476
Murine norovirus, strain S99 Berlin, FLI-RVB-0651	Virucidal - EN 14476

1. Tuffie 5 Disinfecting Wipes

- Bucket of 225
EBOS Code: 10020677
Supplier Code: 901SW225
Bracket
EBOS Code: 10009577
Supplier Code: TTWB225/2275

2. Tuffie 5 Wipes Flexican - 150

EBOS Code: 10020675
Supplier Code: 901SW150FC
Bracket/Holder
EBOS Code: 10009578
Supplier Code: TTWB579/E87579

3. Tuffie 5 Wipes Refill - 150

EBOS Code: 10020676
Supplier Code: 901SW150RP

4. 'Dream' Dispenser

EBOS Code: 10020760
Supplier Code: 520ZS001

5. Tuffie 5 Wipes - 300s

EBOS Code: 10020678
Supplier Code: 901SW300


References: 1. Otter, J.A. et al (2016) Transmission of SARS and MERS coronaviruses and influenza virus in healthcare settings: the possible role of dry surface contamination. Journal of Hospital Infection. 92 235e250 2. Tested to EN 14476 Feline Coronavirus
ARTG 346310

Always read the instructions for use. This product may not be right for you. Read the label before purchase.

TO PLACE AN ORDER, PLEASE GO ONLINE AT www.eboshealthcare.com.au
OR CONTACT US ON **1800 269 534**

ESSENTIAL TO HEALTH

**With
Alcohol**



The **IMPROVED** triple layered fabric, with emboss pockets and microfibre core, absorbs and traps surface soils and holds onto bacteria and bio-film to prevent cross-contamination.

The triple layered construction gives: strength; softness; absorbency; and is non-linting.

The fabric is 100% Polypropylene, optimising bactericidal efficacy.

Tuffie Disinfectant Wipes

Ideal for non porous surfaces where clinical disinfection is required. It contains 70% solution of Isopropyl alcohol BP.

Use in Theatres, A&E, Sterile Services and other areas where required.

Ideal for use on hard surfaces and equipment:




Can be used in Theatres, A&E, Occupational Therapy, Radiology, Catering and anywhere else where clinical cleanliness is required quickly and conveniently.

Other benefits when using Tuffie Disinfectant Wipes:

- 70% solution of Isopropyl Alcohol BP.
- Proven against E.Coli, MRSA plus a range of other micro-organisms.



E.Coli



MRSA

Tuffie Disinfectant Wipes - Canister of 200
EBOS Code: 10020649
Supplier Code: 81106-100

Tuffie Disinfectant Wipes - 6 x 200s
EBOS Code: 10020668
Supplier Code: 901CR200

Bracket
EBOS Code: 10009576
Supplier Code: TTWB100/2515

EBOS HEALTHCARE:
Building 2B, 300 Victoria Street
Wetherill Park, NSW 2164
Free Phone: 1800 269 534
Free Fax: 1800 810 257
Email: customerservice@ebosgroup.com.au
Online Ordering: www.eboshealthcare.com.au

WAREHOUSE LOCATIONS:

New South Wales:
Building 2B, 300 Victoria St
Wetherill Park
NSW 2164

Queensland:
Unit 2, 39 Distribution Street
Motorway Business Park
Larapinta
QLD 4110

Victoria:
34-78 South Park Drive,
South Dandenong,
VIC 3175

Western Australia:
Unit 2, 85 Sheffield Rd
Welshpool
WA 6106

Tasmania:
7A Pear Avenue
Derwent Park
TAS 7009

Vernacare

EBOS
HEALTHCARE

ESSENTIAL TO HEALTH