



WHAT IS VERNAGEL?

Vernagel is designed specifically to cope with bodily fluids, trapping them in a dry semi-solid gel within a bedpan or urinal prior to disposal.

Using Vernagel saves patients the embarrassment and indignity which can follow an accidental spill and reduces the risk of cross-infection to patients and staff.

VERNAGEL

SUPER ABSORBENT POWDER



**SUPER-FAST
ABSORPTION**

- Every sachet contains 7g that can absorb up to 1.2 Litres in total resulting in less spillages, easier disposal and improved infection control.
- Packaging with Vernacare branding and warning message to avoid product mis-use on each individual sachet and dispenser.

Vernacare

FOR MORE INFORMATION CONTACT YOUR ACCOUNT MANAGER
ON **FREEPHONE 1800 269 534** OR VISIT www.eboshealthcare.com.au

EBOS
HEALTHCARE

ESSENTIAL TO HEALTH

WHY VERNAGEL?

- **Reduces the risk of cross infection.** Vernagel solidifies bodily fluids, preventing any spillages when a bedpan or bottle for example is taken from the patient. Cross infection can occur when bodily fluids spill onto beds, floors and other surfaces or equally when receptacles are taken to the dirty utility room for their disposal.
- **Patient dignity.** Vernagel solidifies liquids eliminating spillages and splash backs resulting in increased patient confidence and protecting the dignity of the patient throughout their stay.
- **Healthcare worker satisfaction.** Vernagel eliminates the need to continuously change bedsheets resulting in night shift staff managing their time better. In addition, the minimisation of splashes eradicates any pathogens from being transferred to caregivers when landing on their eyes, nose or mouth. These type of exposures can lead to the acquisition of bloodborne pathogens such as HIV, hepatitis B and C. Occupational exposures such as splashes can result in lost work days, financial burden and can cause emotional distress.
- **Minimises and eliminates residual odours.** Thanks to the solidification process, the potential bad odours are minimised and result in a more pleasant environment for both the healthcare worker and patient.



HOW DOES VERNAGEL WORK?

1. Put **one sachet** of Vernagel into: Bedpans, slipper pans, male or female urinals, vomit bowls, washbowls before use to solidify liquid and prevent spills.
2. **Do not open the sachet.**
3. **Disposal:** Dispose of single-use item and Vernagel into the Vernacare disposal unit. Where not possible, empty into toilet or slop hopper.
Do not tip 'gelled' substance into wash basins or sinks.



**CAN ABSORB
UP TO 1.2 LITRES**



WHERE CAN YOU USE VERNAGEL?

Anywhere where fluids are present:

- In conjunction with Vernacare's single-use containers such as the male and female urinals, bedpans, slipper pans, vomit bowls and washbowls.
- Bins and clinical bags.



**SUPER-FAST
ABSORPTION**

Vernagel large tub

Loose powder with scoop
Weight: 4kg
Supplier Code: 456MA400
EBOS Code: 10020734



Vernagel Bottle

Loose powder
Weight: 475gm
Supplier code: 456MA475
EBOS Code: 10020735



Vernagel 100 Sachet

Easy to use sachets
Weight: 450gm
Supplier Code: 450SC100/450MA100
EBOS Code: 10020733



Vernacare

FOR MORE INFORMATION CONTACT YOUR ACCOUNT MANAGER
ON **FREEPHONE 1800 269 534** OR VISIT www.eboshealthcare.com.au

EBOS
HEALTHCARE

ESSENTIAL TO HEALTH