











**Surgical
Gut – Chronic**

Chromic Gut Suture

Product Characteristics	Chromic Gut Suture
Physical Material	Collagen derived from sheep mucosa and/or beef serosa, treated with chromic salts to slow absorption
Construction	Virtual monofilament
Absorbable / Non-absorbable	Absorbable
Coating (if applicable)	N/A
Colours	Natural – tan (dyed blue for Ophthalmic surgery)
Available Size Range	3 through 7/0
Tensile Strength (Wound Support)	Short to Medium Term
Tensile Strength Retention in Tissue	21-28 days
Absorption Time	Absorbed in 90 days
Frequent Uses	Used in tissues that do not need extended wound support and Ophthalmology

CHROMIC GUT

Surgical Gut Suture – Chromic














▼ 3/8 CIRCLE PRECISION POINT REVERSE CUTTING				5/0 1.5	4/0 2	3/0 3	2/0 3.5	0 4	1 5	2 6
	P-3	13mm	45cm	chromic	687G					
▼ 3/8 CIRCLE REVERSE CUTTING				5/0 1.5	4/0 2	3/0 3	2/0 3.5	0 4	1 5	2 6
	M-2	13mm	30cm	chromic	752G*					
	FS-2	19mm	70cm	chromic	635H	636H				
	FS-1	24mm	70cm	chromic			637H*			
▼ 1/2 CIRCLE PRECISION POINT REVERSE CUTTING				5/0 1.5	4/0 2	3/0 3	2/0 3.5	0 4	1 5	2 6
	PS-5	13mm	45cm	chromic	1653G*					
▼ 1/2 CIRCLE REVERSE CUTTING				5/0 1.5	4/0 2	3/0 3	2/0 3.5	0 4	1 5	2 6
	J-1	17mm	70cm	chromic	U213H^					
	CP-2	26mm	70cm	chromic			887H			
	CP	40mm	70cm	chromic					805H^	
▼ STRAIGHT CUTTING				5/0 1.5	4/0 2	3/0 3	2/0 3.5	0 4	1 5	2 6
	KS	60mm	70cm	chromic			654H*			
▲ 1/2 CIRCLE CONVENTIONAL CUTTING STERNUM				3/0 3	2/0 3.5	0 4	1 5	2 6	3 7	4 8
	CCS-1	60mm	75cm	chromic					W560^	

* New Zealand Indent

^ Australia & New Zealand Indent

CHROMIC GUT

Surgical Gut Suture – Chromic









Needle length and description		Length	Colour	Gauge						
⊙ 3/8 CIRCLE TAPER POINT				4/0 2	3/0 3	2/0 3.5	0 4	1 5	2 6	3 7
BB 	17mm	70cm	chromic	4887G						
⊙ 1/2 CIRCLE TAPER POINT				4/0 2	3/0 3	2/0 3.5	0 4	1 5	2 6	3 7
RB-1 	17mm	70cm	chromic	U203H	U204H*					
CT-3 	22mm	70cm	chromic		892H^	893H^				
SH-1 	22mm	70cm	chromic	G181H	G182H					
SH 	26mm	70cm	chromic	G121H	G122H	G123H	G124H*			
SH 	26mm 26mm	90cm	chromic			G173H^				
CT-2 	27mm	70cm	chromic			883H^	884H*	885H^		
MO-5 	31mm	70cm	chromic					4995H*		
CT-1 	36mm	70cm 90cm	chromic chromic		810H^	811H*	812H	813H 925H	814H	
MH 	36mm	70cm	chromic		G126H^	G127H	G128H			
MO-4 	36mm	70cm	chromic				4924H*	4925H		
CT 	40mm	70cm 90cm	chromic chromic			801H 913H^	802H* 914H^	803H 915H*		
LH 	40mm	70cm	chromic			G133H^				

* New Zealand Indent

^ Australia & New Zealand Indent

CHROMIC GUT

Surgical Gut Suture – Chromic






Needle length and description		Length	Colour	Gauge						
				4/0 2	3/0 3	2/0 3.5	0 4	1 5	2 6	3 7
☉ 1/2 CIRCLE TAPER POINT										
 MO-2	40mm	70cm	chromic					4915H		
 CTX	48mm	70cm 90cm	chromic chromic			863H*	904H*	865H 905H		
☽ 1/2 CIRCLE TAPER POINT, 1/2 CIRCLE REVERSE CUTTING										
 CT-1 / CP-1	36mm 36mm	150cm	chromic			983H*				
☉ 5/8 CIRCLE TAPER POINT										
 UR-6	27mm	70cm	chromic			N878H^				
☾ 1/2 CIRCLE TAPERCUT										
 V-7	26mm	75cm	chromic			W565^				
 V-34	36mm	90cm	chromic			943H^	944H^	945H^		
▲ TROCAR POINT										
 V-38	31mm	75cm	chromic					W928^		
 V-37	40mm	75cm	chromic					W940*		

* New Zealand Indent

^ Australia & New Zealand Indent

CHROMIC GUT

Surgical Gut Suture – Chromic

Needle length and description	Length	Colour	Gauge						
▲ TROCAR POINT			4/0 2	3/0 3	2/0 3.5	0 4	1 5	2 6	3 7
 V-40 48mm	75cm	chromic							W947 [^]
▼ 3/8 CIRCLE MICROPOINT REVERSE CUTTING			6/0 1	5/0 1.5	4/0 2	3/0 3	2/0 3.5	0 4	1 5
 G-1 10.5mm	45cm	chromic	796G						
 G-1 10.5mm 10.5mm	45cm	chromic	790G*						
▾ 1/4 CIRCLE SABRELOC CONVENTIONAL SPATULA			7/0 0.5	6/0 1	5/0 1.5	4/0 2	3/0 3	2/0 3.5	0 4
 S-4 8mm 8mm	45cm	chromic				1793G*			
▾ 1/4 CIRCLE SABRELOC SPATULA			7/0 0.5	6/0 1	5/0 1.5	4/0 2	3/0 3	2/0 3.5	0 4
 S-14 8mm 8mm	45cm 45cm	chromic blue		1791G*	1792G*				
NON-NEEDED			4/0 2	3/0 3	2/0 3.5	0 4	1 5	2 6	3 7
LIGAPAK*	135cm	chromic			L113G	L114G*	L115G*		
NON-NEEDED	150cm	chromic	S111H[^]	S112H*	S113H	S114H	S115H*	S116H*	

* New Zealand Indent

[^] Australia & New Zealand Indent