

The background is a solid teal color. On the left side, there are several white, semi-transparent, flowing lines that resemble smoke or liquid trails, moving downwards and slightly to the right. The word "ETHILON" is written in a bold, white, sans-serif font, positioned in the middle-right area of the image.




ETHILON®








ETHILON* Nylon Suture




Product Characteristics	ETHILON* Suture
Physical Material	Polyamide polymer and Nylon 6
Construction	Monofilament
Absorbable / Non-absorbable	Non-absorbable
Coating (if applicable)	N/A
Colours	Clear, black and green
Available Size Range	2 through 11/0
Tensile Strength (Wound Support)	Long term wound support
Tensile Strength Retention in Tissue	Non-absorbable, but loses about 10-15% strength per year
Absorption Time	Degrades at a rate of about 10-15% per year
Frequent Uses	Skin closure; retention; Plastic surgery; Ophthalmology and Microsurgery

*Trademark

ETHILON Nylon

Needle length and description		Length	Colour	Gauge						
○ 1/2 CIRCLE ETHIGUARD BLUNT POINT				0 3.5	1 4	2 5	3 6	4 6	5 7	7 9
	36mm	100cm	black		B020G					
	48mm	150cm loop	black		W729					
	65mm	240cm loop	black		D8555^					





▼ 3/8 CIRCLE PRECISION POINT REVERSE CUTTING				7/0 0.5	6/0 0.7	5/0 1	4/0 1.5	3/0 2	2/0 3	0 3.5
	11mm	45cm 45cm	black black undyed	1696G	697G W1610T^ 689G^					
	13mm	45cm 45cm 45cm	black black undyed			698G W1611T^ 690G*	699G^ 691G*			
	13mm	45cm	black		1698G					
	16mm	45cm	black		1665G*	1668G				
	19mm	45cm 75cm	undyed undyed				1611G^		D9635	
	19mm	45cm	black			1666G	1667G			
	24mm	45cm 75cm	black undyed					1663H D9782^		



▼ 3/8 CIRCLE REVERSE CUTTING				7/0 0.5	6/0 0.7	5/0 1	4/0 1.5	3/0 2	2/0 3	0 3.5
	13mm	45cm	black		667G	668G				
	19mm	45cm	black				662G	653G		
	19mm	45cm	black			661G				


* New Zealand Indent


^ Australia & New Zealand Indent

ETHILON Nylon

Needle length and description		Length	Colour	Gauge						
▼ 3/8 CIRCLE REVERSE CUTTING				7/0 0.5	6/0 0.7	5/0 1	4/0 1.5	3/0 2	2/0 3	0 3.5
 FS-1	24mm	45cm 75cm	black black				1629H	663G 669H		
 FS	26mm	45cm	black						664G	
 FSL	30mm	75cm	black				1670H*	1671G		
 FSLX	40mm	75cm	black					1673H	1674H	

▼ 3/8 CIRCLE REVERSE CUTTING				0 3.5	1 4	2 5	3 6	4 6	5 7	7 9
 LR	75mm	50cm	black			460T^				
 LR	75mm 75mm	75cm	black			489T^				

▼ 1/2 CIRCLE PRECISION POINT REVERSE CUTTING				6/0 0.7	5/0 1	4/0 1.5	3/0 2	2/0 3	0 3.5	1 4
 PS-4	16mm	45cm 45cm	black undyed			1662G 1603G				



▼ 1/2 CIRCLE REVERSE CUTTING				6/0 0.7	5/0 1	4/0 1.5	3/0 2	2/0 3	0 3.5	1 4
 X-1	22mm	45cm	black				642G*			


PLIABILISED







* New Zealand Indent


^ Australia & New Zealand Indent

ETHILON Nylon

Needle length and description		Length	Colour	Gauge						
▼ 1/2 CIRCLE REVERSE CUTTING				6/0 0.7	5/0 1	4/0 1.5	3/0 2	2/0 3	0 3.5	1 4
	40mm	100cm 150cm loop	black black						W737 D10109^	
	48mm	100cm	black							W738*

▼ STRAIGHT CUTTING				6/0 0.7	5/0 1	4/0 1.5	3/0 2	2/0 3	0 3.5	1 4
	60mm	75cm	black			626H*	627H	628H		

▲ 3/8 CIRCLE CONVENTIONAL CUTTING PRIME				6/0 0.7	5/0 1	4/0 1.5	3/0 2	2/0 3	0 3.5	1 4
	13mm	45cm	black	1856G						
	13mm	45cm	black	1956G*	1955G*					
PLIABILISED										
	16mm	45cm	black	W1615T	1865G					
	16mm	45cm	black	1966G*	1965G*	1964G				
PLIABILISED										
	19mm	45cm	black			1894G	1893G			
	19mm	45cm	black			1994G*	1993G*			
PLIABILISED										





▲ 3/8 CIRCLE CONVENTIONAL CUTTING				6/0 0.7	5/0 1	4/0 1.5	3/0 2	2/0 3	0 3.5	1 4
	16mm	35cm	black	W506*						





* New Zealand Indent

^ Australia & New Zealand Indent

ETHILON Nylon

Needle length and description	Length	Colour	Gauge						
▲ 3/8 CIRCLE CONVENTIONAL CUTTING			2/0 3	0 3.5	1 4	2 5	3 6	4 6	5 7
 <p style="text-align: center;">CKP-3</p> <p style="text-align: right;">90mm</p> <p>WITH 40MM RETENTION TUBING</p>	2x50cm	black			W2797				











◎ 3/8 CIRCLE TAPER POINT				11/0 0.1	10/0 0.2	9/0 0.3	8/0 0.4	7/0 0.5	6/0 0.7	5/0 1
BV130-3 	3.8mm	12cm	black		2820G*	2819G*				
BV75-3 	3.8mm	12cm	black		2870G					
BV50-3 	4.1mm	13cm 13cm	black black	2881G* W2881^						
BV130-4 	5mm	12cm 15cm	black black		2814G*	2813G	2815G			
BV100-4 	5.1mm	15cm	black		2830G^	2829G^				
BV75-4 	5.1mm	12cm	black		2850G*					
BV130-5 	6.5mm	12cm 15cm	black black		2810G*	2809G	2808G			
BV-1 	9.3mm 9.3mm	45cm	black				W8170^			

◎ 1/2 CIRCLE TAPER POINT				0 3.5	1 4	2 5	3 6	4 6	5 7	7 9
MH-1 	31mm	100cm	black	4784G*	W745					
MH-1 	31mm	150cm loop	black	W739^						
CT-1 	36mm	150cm loop	black		W743*					
CT 	40mm	100cm 100cm	black black	W746^	4999G W747					

* New Zealand Indent

^ Australia & New Zealand Indent























ETHILON Nylon

Needle length and description		Length	Colour	Gauge						
☉ 1/2 CIRCLE TAPER POINT				0 3.5	1 4	2 5	3 6	4 6	5 7	7 9
	40mm	150cm loop	black	W740						
	40mm	200cm loop	black	W760^						
	48mm	100cm	black		4985G	4986G*				
	48mm	150cm loop	black	W741*	W748					
	48mm	200cm loop	black		W768^					
	65mm	120cm loop	black	L880G^						
	65mm	150cm	black			825G^				
☽ 1/4 CIRCLE TAPERCUT				10/0 0.2	9/0 0.3	8/0 0.4	7/0 0.5	6/0 0.7	5/0 1	4/0 1.5
	3.5mm	15cm	black		2894G^					
	3.5mm	15cm	black	2892G*						
	3.8mm	15cm	black			2897G*				

* New Zealand Indent

^ Australia & New Zealand Indent









ETHILON Nylon

Needle length and description		Length	Colour	Gauge						
 1/4 CIRCLE MICROPOINT SPATULA				10/0 0.2	9/0 0.3	8/0 0.4	7/0 0.5	6/0 0.7	5/0 1	4/0 1.5
 TG100-8 6.5mm 6.5mm	30cm	black			1714G*					
 3/8 CIRCLE MICROPOINT SPATULA				11/0 0.1	10/0 0.2	9/0 0.3	8/0 0.4	7/0 0.5	6/0 0.7	5/0 1
 VAS100-4 5.1mm	15cm	black			2890G*					
 TG140-6 6.5mm 6.5mm	30cm	black	W1780^	W1770^	W1769^					
 TG140-8 6.5mm 6.5mm	30cm	black		7718G						
 1/2 CIRCLE MICROPOINT SPATULA				11/0 0.1	10/0 0.2	9/0 0.3	8/0 0.4	7/0 0.5	6/0 0.7	5/0 1
 TG160-6 5.5mm 5.5mm	15cm 30cm	black black		7757G^	W1756*					
 TG175-8 7mm 7mm	30cm	black				1716G*				
 1/4 CIRCLE SABRELOC SPATULA				11/0 0.1	10/0 0.2	9/0 0.3	8/0 0.4	7/0 0.5	6/0 0.7	5/0 1
 S-24 8mm 8mm	45cm	black								7731G
 1/4 CIRCLE SABRELOC CENTRE POINT SPATULA				11/0 0.1	10/0 0.2	9/0 0.3	8/0 0.4	7/0 0.5	6/0 0.7	5/0 1
 RD-1 8mm 8mm	45cm	black								749G
<small>PLIABILISED</small>										
 1/4 CIRCLE ULTIMA SPATULA				11/0 0.1	10/0 0.2	9/0 0.3	8/0 0.4	7/0 0.5	6/0 0.7	5/0 1
 CS100-8 6.7mm 6.7mm	30cm	black				9039G				
 CS90-6 7mm	10cm	black		9800G^						
 3/8 CIRCLE ULTIMA SPATULA				11/0 0.1	10/0 0.2	9/0 0.3	8/0 0.4	7/0 0.5	6/0 0.7	5/0 1
 CS140-6 6.5mm	10cm	black		9060G						
 CS140-6 6.5mm 6.5mm	30cm	black		9003G						
 CS140-8 6.75mm 6.75mm	30cm	black		9033G^						
 1/2 CIRCLE ULTIMA SPATULA				11/0 0.1	10/0 0.2	9/0 0.3	8/0 0.4	7/0 0.5	6/0 0.7	5/0 1
 CSM-6 5mm	30cm	black		9005G^						

* New Zealand Indent

^ Australia & New Zealand Indent

ETHILON Nylon

Needle length and description		Length	Colour	Gauge						
 1/2 CIRCLE ULTIMA SPATULA				11/0 0.1	10/0 0.2	9/0 0.3	8/0 0.4	7/0 0.5	6/0 0.7	5/0 1
CS160-6 	5.5mm	10cm	black		9061G					
CS160-6 	5.5mm 5.5mm	15cm	black		9001G*					
CS160-6 	5.5mm	20cm	black	9016G^						
CS160-6 	5.5mm 5.5mm	30cm	black		9000G*					
CS160-8 	5.5mm 5.5mm	30cm	black		9030G^					
CS175-6 	5.5mm 5.5mm	30cm	black		9002G^					
CS175-8 	7mm 7mm	30cm	black		9032G*	9036G	9038G^			
NON-NEEDED				3/0 2	2/0 3	0 3.5	1 4	2 5	3 6	4 6
NON-NEEDED		2x100cm	black				W185			

* New Zealand Indent

^ Australia & New Zealand Indent

