



MERSILENE®


MERSILENE* Polyester Suture

Product Characteristics	MERSILENE* Suture
Physical Material	Polyester fibre
Construction	Braided (Monofilament sizes 10/0 and 11/0)
Absorbable / Non-absorbable	Non-absorbable
Coating (if applicable)	N/A
Colours	White and green
Available Size Range	5 through 11/0
Tensile Strength (Wound Support)	Long term wound support
Tensile Strength Retention in Tissue	Indefinite
Absorption Time	Non-absorbable
Frequent Uses	Cardiovascular; General and Plastic surgery; retention; Ophthalmology

*Trademark

MERSILENE

Polyester Fiber Suture



Needle length and description	Length	Colour	Gauge						
○ 1/2 CIRCLE ETHIGUARD BLUNT POINT			5mm tape						
 BP-1 65mm 65mm	30cm	undyed	RS21						
▼ 3/8 CIRCLE PRECISION POINT REVERSE CUTTING			5/0 1	4/0 1.5	3/0 2	2/0 3	0 3.5	1 4	2 5
 P-3 13mm	45cm	undyed	R690G^	R691G					
▼ 3/8 CIRCLE REVERSE CUTTING			5/0 1	4/0 1.5	3/0 2	2/0 3	0 3.5	1 4	2 5
 FS-2 19mm	45cm	undyed		R633H*					
 FS-1 24mm	45cm	undyed			R647H*				
▼ 1/2 CIRCLE REVERSE CUTTING			5/0 1	4/0 1.5	3/0 2	2/0 3	0 3.5	1 4	2 5
 CP 40mm	75cm	undyed						R995G	
▼ 1/2 CIRCLE REVERSE CUTTING			5mm tape						
 OS-8 40mm	90cm	undyed	D9212^						
⊙ 1/2 CIRCLE TAPER POINT			3/0 2	2/0 3	0 3.5	1 4	2 5	3 6	4 6
 SH 26mm	75cm	green	R832H^	R833H*					

* New Zealand Indent

^ Australia & New Zealand Indent

MERSILENE

Polyester Fiber Suture

Needle length and description	Length	Colour	Gauge						
 1/2 CIRCLE TAPER POINT			5mm tape						
 CTX 48mm 48mm	30cm	undyed	RS22						
 STRAIGHT CUTTING			5/0 1	4/0 1.5	3/0 2	2/0 3	0 3.5	1 4	2 5
 KS 60mm	75cm	white				D9433^			
 1/2 CIRCLE MICROPOINT SPATULA			11/0 0.1	10/0 0.2	9/0 0.3	8/0 0.4	7/0 0.5	6/0 0.7	5/0 1
 TG160-4 5.5mm 5.5mm	30cm	green	R750G^						
 1/4 CIRCLE SABRELOC CONVENTIONAL SPATULA			4/0 1.5	3/0 2	2/0 3	0 3.5	1 4	2 5	3 6
 S-4 8mm 8mm	45cm	undyed	D2220						
 1/2 CIRCLE SABRELOC CONVENTIONAL SPATULA			6/0 0.7	5/0 1	4/0 1.5	3/0 2	2/0 3	0 3.5	1 4
 S-28 8mm	45cm	undyed	D7683^						
 1/4 CIRCLE SABRELOC SPATULA			6/0 0.7	5/0 1	4/0 1.5	3/0 2	2/0 3	0 3.5	1 4
 S-14 8mm 8mm	45cm	undyed		1760G					
 S-14 8mm 8mm	2x45cm	undyed		W843	W832*				
 1/2 CIRCLE SABRELOC SPATULA			6/0 0.7	5/0 1	4/0 1.5	3/0 2	2/0 3	0 3.5	1 4
 S-2 8mm 8mm	45cm	undyed			1779G				
 1/4 CIRCLE SABRELOC CENTRE POINT SPATULA			6/0 0.7	5/0 1	4/0 1.5	3/0 2	2/0 3	0 3.5	1 4
 RD-1 8mm 8mm	45cm	undyed		746G^					
 3/8 CIRCLE ULTIMA SPATULA			10/0 0.2	9/0 0.3	8/0 0.4	7/0 0.5	6/0 0.7	5/0 1	4/0 1.5
 CS140-6 6.5mm 6.5mm	30cm	green	R770G^						

* New Zealand Indent

^ Australia & New Zealand Indent

MERSILENE

Polyester Fiber Suture

Needle length and description		Length	Colour	Gauge						
 1/2 CIRCLE ULTIMA SPATULA				10/0 0.2	9/0 0.3	8/0 0.4	7/0 0.5	6/0 0.7	5/0 1	4/0 1.5
CS160-6 	5.5mm 5.5mm	30cm	green	R756G^						

 BI-CURVE ULTIMA SPATULA				10/0 0.2	9/0 0.3	8/0 0.4	7/0 0.5	6/0 0.7	5/0 1	4/0 1.5
CSB-6 	4.8mm 4.8mm	20cm	green	R721G^						

