

















MONOCRYL®

MONOCRYL* (poliglecaprone 25) Suture

Product Characteristics	MONOCRYL* Suture
Physical Material	Poliglecaprone 25 (copolymer of glycolide and epsilon-caprolactone)
Construction	Monofilament
Absorbable / Non-absorbable	Absorbable
Coating (if applicable)	N/A
Colours	Undyed – beige and Dyed – violet
Available Size Range	2 through 6/0
Tensile Strength (Wound Support)	Medium Term (3-4 weeks)
Tensile Strength Retention in Tissue	Undyed - 50-60% of original strength remains at 1 week - 20-30% of original strength remains at 2 weeks Dyed - 60-70% of original strength remains at 1 week - 30-40% of original strength remains at 2 weeks
Absorption Time	Completely absorbed in 91-119 days

*Trademark




MONOCRYL® (Poliglecaprone 25) Monofilament Suture











Needle length and description		Length	Colour	Gauge						
○ 1/2 CIRCLE ETHIGUARD BLUNT POINT				3/0 2	2/0 3	0 3.5	1 4	2 5	3 6	4 6
 SHB	26mm	70cm	undyed	YB416*						
 CTB-1	36mm	70cm 90cm	violet undyed	YB944^		YB258* YB946^	YB947*			
 CTX-B	48mm	90cm	undyed			YB978^				
 XLH-D	70mm	90cm	violet				W3709			
▼ 3/8 CIRCLE PRECISION POINT REVERSE CUTTING				6/0 0.7	5/0 1	4/0 1.5	3/0 2	2/0 3	0 3.5	1 4
 P-1	10.5mm	45cm	undyed	Y489G	Y490G					
 P-3	13mm	45cm	undyed	Y492G	Y493G	Y494G				
 PS-2	19mm	45cm 70cm	undyed undyed			Y496G Y426H	Y427H			
 PS-1	24mm	70cm 90cm	undyed undyed			Y935H D8797	Y936H			
 PS	26mm	70cm	undyed				W3213			
▼ 3/8 CIRCLE REVERSE CUTTING				6/0 0.7	5/0 1	4/0 1.5	3/0 2	2/0 3	0 3.5	1 4
 FS	26mm	70cm	undyed				W3327			
▼ 1/2 CIRCLE PRECISION POINT REVERSE CUTTING				6/0 0.7	5/0 1	4/0 1.5	3/0 2	2/0 3	0 3.5	1 4
 P-2	8mm	45cm	undyed			D8539^				
 PS-6	10.5mm	45cm	undyed	Y510G						
▼ 1/2 CIRCLE REVERSE CUTTING				6/0 0.7	5/0 1	4/0 1.5	3/0 2	2/0 3	0 3.5	1 4
 CP-1	36mm	70cm	undyed						Y267H	
▼ STRAIGHT CUTTING				6/0 0.7	5/0 1	4/0 1.5	3/0 2	2/0 3	0 3.5	1 4
 KS	60mm	70cm	undyed				Y523H			

* New Zealand Indent

^ Australia & New Zealand Indent

MONOCRYL® (Poliglecaprone 25) Monofilament Suture





Needle length and description		Length	Colour	Gauge						
▲ 3/8 CIRCLE CONVENTIONAL CUTTING PRIME				6/0 0.7	5/0 1	4/0 1.5	3/0 2	2/0 3	0 3.5	1 4
 PC-1	13mm	45cm	undyed	Y833G	Y834G	Y835G				
 PC-3	16mm	45cm	undyed		Y844G	Y845G				
 PC-5	19mm	45cm	undyed		Y822G*	Y823G				

◎ 1/2 CIRCLE TAPER POINT				6/0 0.7	5/0 1	4/0 1.5	3/0 2	2/0 3	0 3.5	1 4
 TF	13mm	70cm	undyed	Y432H						
 RB-1	17mm	70cm 70cm	violet undyed		Y303H	Y304H	Y305H Y215H*			
 JB	26mm	70cm	violet					W3665^		
 SH	26mm	70cm 70cm	violet undyed			Y315H	Y316H	Y317H Y417H*		
 SH	26mm 26mm	90cm	violet				Y527H			
 CT-2	27mm	70cm	violet					Y333H	Y334H	Y335H*
 CT-1	36mm	90cm	violet					Y345H	Y346H	Y347H
 MH	36mm	70cm 70cm	undyed violet					Y227H* W3463		
 CT	40mm	90cm	violet					Y357H*	Y358H*	Y359H
 CTX	48mm	90cm	violet					Y397H^		Y399H*

* New Zealand Indent

^ Australia & New Zealand Indent

MONOCRYL® (Poliglecaprone 25) Monofilament Suture

Needle length and description	Length	Colour	Gauge						
⊙ 5/8 CIRCLE TAPER POINT			6/0 0.7	5/0 1	4/0 1.5	3/0 2	2/0 3	0 3.5	1 4
 UR-6 27mm	70cm	violet				Y604H	Y605H	Y606H*	
⊙ 5/8 CIRCLE TAPER POINT			6/0 0.7	5/0 1	4/0 1.5	3/0 2	2/0 3	0 3.5	1 4
 UR-6 27mm 27mm	70cm 90cm	undyed violet				D9645 D9503	D8592		
▽ 1/2 CIRCLE TAPERCUT			6/0 0.7	5/0 1	4/0 1.5	3/0 2	2/0 3	0 3.5	1 4
 V-34 36mm	90cm	undyed					Y917H*		
⊙ 1/2 CIRCLE TAPER POINT, 1/2 CIRCLE REVERSE CUTTING			6/0 0.7	5/0 1	4/0 1.5	3/0 2	2/0 3	0 3.5	1 4
 CT-1 / CP-1 36mm 36mm	135cm	undyed				Y931H^	Y932H^		
NON-NEEDED			6/0 0.7	5/0 1	4/0 1.5	3/0 2	2/0 3	0 3.5	1 4
NON-NEEDED	6x45cm	violet					Y909G*	Y910G^	Y911G*

* New Zealand Indent

^ Australia & New Zealand Indent

