

The background is a solid teal color. Overlaid on this are several thin, white, semi-transparent lines that flow and curve across the page, creating a sense of movement and depth. These lines are most prominent on the left side and extend towards the center.

NUROLON®










NUROLON* Braided Nylon Suture

Product Characteristics	NUROLON* Suture
Physical Material	Nylon 6
Construction	Braided
Absorbable / Non-absorbable	Non-absorbable
Coating (if applicable)	N/A
Colours	Black and white
Available Size Range	1 through 6/0
Tensile Strength (Wound Support)	Long term wound support
Tensile Strength Retention in Tissue	Non-absorbable, but loses about 10-15% strength per year
Absorption Time	Degrades at a rate of about 10-15% per year
Frequent Uses	Most body tissues for ligating and suturing; General closure; Neurosurgery

*Trademark

NUROLON

Braided Nylon Suture

Needle length and description		Length	Colour	Gauge						
▼ 1/2 CIRCLE REVERSE CUTTING				3/0 2	2/0 3	0 3.5	1 4	2 5	3 6	4 6
 X-1	22mm	8x45cm control release	black	C563D [^]						
▼ STRAIGHT CUTTING				3/0 2	2/0 3	0 3.5	1 4	2 5	3 6	4 6
 KS	60mm	75cm	black	5622H [^]						
▲ 3/8 CIRCLE CONVENTIONAL CUTTING PRIME				5/0 1	4/0 1.5	3/0 2	2/0 3	0 3.5	1 4	2 5
 PC-5	19mm	45cm	black		W5511 [^]					
⊙ 1/2 CIRCLE TAPER POINT				5/0 1	4/0 1.5	3/0 2	2/0 3	0 3.5	1 4	2 5
 TF	13mm	8x45cm control release	black		C584D [*]					
 TF	13mm	8x45cm	black		N104T [*]					
 RB-1	17mm	8x45cm control release	black		C554D [^]	C553D [*]				
 RB-1	17mm	8x45cm	black		N124T [*]					
 RB-1	17mm	75cm	black		D1637					
 SH-1	22mm	8x45cm control release	black			C503D [*]				

* New Zealand Indent

[^] Australia & New Zealand Indent

