



PDS™


PDS™ II (polydioxanone) Suture

Product Characteristics	PDS™ II Suture
Physical Material	Polydioxanone
Construction	Monofilament
Absorbable / Non-absorbable	Absorbable
Coating (if applicable)	N/A
Colours	Undyed – clear and Dyed – violet
Available Size Range	2 through 9/0
Tensile Strength (Wound Support)	Long Term (8 weeks)
Tensile Strength Retention in Tissue	4/0 and smaller: 60% of original strength remains at 2 weeks 40% of original strength remains at 4 weeks 35% of original strength remains at 6 weeks 3/0 and larger: 80% of original strength remains at 2 weeks 70% of original strength remains at 4 weeks 60% of original strength remains at 6 weeks
Absorption Time	Completely absorbed in 182-238 days
Frequent Uses	Abdominal and thoracic closure; subcutaneous tissue; Colon and Rectal surgery; Orthopaedic and Plastic Surgery

*Trademark

PDS™ II

(Polydioxanone) Monofilament Suture

Needle length and description		Length	Colour	Gauge						
○ 1/2 CIRCLE ETHIGUARD BLUNT POINT				0 3.5	1 4	2 5	3 6	4 6	5 7	7 9
 CTB-1	36mm	70cm 90cm	violet violet	ZB340	ZB347*					
 TD-3	45mm	90cm 90cm	violet violet		W9965T					
 CTX-B	48mm	90cm	violet	ZB370^	ZB371					
 CTX-D	48mm	150cm loop	violet		W9967T					
 BP-1	65mm	244cm loop	violet		D8892					

▼ 3/8 CIRCLE PRECISION POINT REVERSE CUTTING				6/0 0.7	5/0 1	4/0 1.5	3/0 2	2/0 3	0 3.5	1 4
 P-1	11mm	45cm	undyed	Z489G						
 P-3	13mm	45cm 45cm	violet undyed		Z463G* Z493G	Z494G^				
 PS-3	16mm	45cm	undyed			W9873T^				
 PS-2	19mm	45cm 45cm 45cm	violet undyed undyed			Z513G* Z496G W9874T	Z497G			
 PS-1	24mm	45cm	undyed				Z683G			

▼ 3/8 CIRCLE REVERSE CUTTING				6/0 0.7	5/0 1	4/0 1.5	3/0 2	2/0 3	0 3.5	1 4
 FS-2	19mm	45cm 70cm	undyed undyed		Z421H*	W9950T^ Z422H	Z423H			

* New Zealand Indent



^ Australia & New Zealand Indent

PDS™ II




(Polydioxanone) Monofilament Suture

Needle length and description	Length	Colour	Gauge							
▼ 3/8 CIRCLE REVERSE CUTTING			6/0 0.7	5/0 1	4/0 1.5	3/0 2	2/0 3	0 3.5	1 4	
 FS-1	24mm	70cm	undyed			Z441H	Z442H	Z443H		
 FS	26mm	45cm	undyed				W9957T^			

▼ 1/2 CIRCLE PRECISION POINT REVERSE CUTTING			6/0 0.7	5/0 1	4/0 1.5	3/0 2	2/0 3	0 3.5	1 4	
 P-2	8mm	45cm	undyed		Z503G^	Z504G^				
 PS-4	16mm	45cm	undyed			Z507G				

▼ 1/2 CIRCLE REVERSE CUTTING			3/0 2	2/0 3	0 3.5	1 4	2 5	3 6	4 6	
 OS-4	22mm	70cm 70cm	violet violet				W9286T*			
 CP-2	26mm	70cm	violet	Z968H^	Z969H	Z970H*				
 CP-1	36mm	70cm	violet		Z466H*	Z467H*	Z468H			
 OS-6	36mm	70cm	violet				Z535H*			
 CP	40mm	70cm	violet					Z195T		

▼ STRAIGHT CUTTING			4/0 1.5	3/0 2	2/0 3	0 3.5	1 4	2 5	3 6	
 KS	60mm	70cm 70cm	violet undyed	Z662H^	Z663H^	W9715T^				

▲ 3/8 CIRCLE CONVENTIONAL CUTTING PRIME			6/0 0.7	5/0 1	4/0 1.5	3/0 2	2/0 3	0 3.5	1 4	
 PC-1	13mm	45cm	undyed	Z833G						
 PC-3	16mm	45cm	undyed		Z844G					
 PC-5	19mm	45cm 45cm	violet undyed			W9615T Z823G	Z824G			
 PC-25	26mm	45cm 70cm	undyed undyed			W9870T	W9625T			

* New Zealand Indent

^ Australia & New Zealand Indent

PDS™ II

(Polydioxanone) Monofilament Suture



3/8 CIRCLE TAPER POINT				7/0 0.5	6/0 0.7	5/0 1	4/0 1.5	3/0 2	2/0 3	0 3.5
BV-1 	9.3mm 9.3mm	60cm 75cm	violet violet	Z155H^ Z135H^	Z117H					
BV 	10.5mm 10.5mm	75cm	violet		Z123H^					
C-1 	13mm 13mm	75cm	violet		Z127H^					


1/2 CIRCLE TAPER POINT				6/0 0.7	5/0 1	4/0 1.5	3/0 2	2/0 3	0 3.5	1 4
RB-2 	13mm 13mm	75cm	violet	Z149H*						
TF 	13mm	70cm	violet	Z432H						
UCL 	14mm	70cm	violet					Z114H^		
RB-1 	17mm	70cm	violet		Z303H	Z304H	Z305H			
RB-1 	17mm 17mm	75cm 90cm	violet violet		Z320H	W9109H*				
SH-1 	22mm	70cm	violet			Z310H^				
SH 	26mm	8x45cm control release	violet			Z773D^	Z774D*	Z775D^		
SH 	26mm	70cm	violet			Z315H	Z316H	Z317H		
SH 	26mm 26mm	90cm	violet				D6264*			
CT-2 	27mm	70cm	violet					Z333H	Z334H	
CT-3 	22mm	70cm	violet						Z329H	

* New Zealand Indent

^ Australia & New Zealand Indent

PDS™ II (Polydioxanone) Monofilament Suture

Needle length and description		Length	Colour	Gauge						
☉ 1/2 CIRCLE TAPER POINT				6/0 0.7	5/0 1	4/0 1.5	3/0 2	2/0 3	0 3.5	1 4
 CT-1	36mm	70cm 90cm	violet violet				Z338H	Z339H	Z340H Z346H*	Z341H Z347H
 CT	40mm	70cm 90cm	violet violet					Z357H^	Z352H* Z358T	Z353H* Z359T
 MO-2	40mm	150cm loop	violet						W9236T^	
 CTX	48mm	90cm	violet						Z370T^	Z371T
 CTX	48mm	150cm loop	violet							W9262T^
 TP-1	65mm	75cm 150cm loop 244cm loop	violet violet violet						Z991G*	D7819 Z880G
 XLH	70mm	70cm 135cm	violet violet							Z583G^ Z877G^

☉ 5/8 CIRCLE TAPER POINT				6/0 0.7	5/0 1	4/0 1.5	3/0 2	2/0 3	0 3.5	1 4
 UR-6	27mm	70cm	violet						D7185	
 MF-110	110mm	150cm loop	violet						W9168^	

* New Zealand Indent

^ Australia & New Zealand Indent

PDS™ II

(Polydioxanone) Monofilament Suture

Needle length and description		Length	Colour	Gauge						
☉ STRAIGHT TAPER POINT				6/0 0.7	5/0 1	4/0 1.5	3/0 2	2/0 3	0 3.5	1 4
 ST-4	19mm 19mm	135cm	violet			Z420G*				

☽ 1/2 CIRCLE TAPERCUT				6/0 0.7	5/0 1	4/0 1.5	3/0 2	2/0 3	0 3.5	1 4
 V-38	31mm	70cm	violet						W9333T*	
 V-34	35mm	90cm	violet						W9381T^	W9352T^
 V-37	40mm	240cm loop	violet							W9394

☾ STRAIGHT TROCAR POINT				6/0 0.7	5/0 1	4/0 1.5	3/0 2	2/0 3	0 3.5	1 4
 STP-10	254mm 254mm	70cm	violet					Z997G*	D9864^	

* New Zealand Indent

^ Australia & New Zealand Indent

