



**PROLENE®**





## PROLENE\* Polypropylene Suture







Product Characteristics	PROLENE* Suture
Physical Material	Polypropylene
Construction	Monofilament
Absorbable / Non-absorbable	Non-absorbable
Coating (if applicable)	N/A
Colours	Clear and blue
Available Size Range	2 through 10/0
Tensile Strength (Wound Support)	Long term wound support
Tensile Strength Retention in Tissue	Indefinite
Absorption Time	Non-absorbable
Frequent Uses	General; Plastic; Cardiovascular surgery; skin closure and Ophthalmology




\*Trademark

# PROLENE

## Monofilament Polypropylene Suture

Needle length and description	Length	Colour	Gauge						
<b>○ 1/2 CIRCLE ETHIGUARD BLUNT POINT</b>			0 3.5	1 4	2 5	3 6	4 6	5 7	7 9
 CTB-1 36mm	75cm 75cm	undyed blue	<b>D9776*</b>						
 CTD 40mm	100cm	blue		<b>W8998*</b>					
 CTX-B 48mm	75cm	blue	<b>PB31H^</b>	<b>PB30H</b>					
 BP-1 65mm	150cm loop	blue		<b>D9668</b>					

▼ 3/8 CIRCLE PRECISION POINT REVERSE CUTTING			7/0 0.5	6/0 0.7	5/0 1	4/0 1.5	3/0 2	2/0 3	0 3.5
 P-1 10.5mm	45cm 45cm	blue undyed	<b>8696G^</b>	<b>8697G^</b> <b>8606G^</b>					
 P-3 13mm	45cm	blue			<b>8698G</b>				
 PS-3 16mm	45cm	blue		<b>8680G^</b>	<b>8681G</b>				
 PS-2 19mm	45cm 45cm	blue blue			<b>8686G</b> <b>W8010T^</b>	<b>8682H</b> <b>W8015T</b>	<b>8687H*</b>		
 PS-1 24mm	45cm	blue					<b>8663G</b>		
 PS 26mm	45cm	blue					<b>W8025T</b>		











▼ 3/8 CIRCLE REVERSE CUTTING			5/0 1	4/0 1.5	3/0 2	2/0 3	0 3.5	1 4	2 5
 FS-3 16mm	45cm	blue	<b>8614H*</b>						
 FS-2 19mm	45cm	blue	<b>8661G*</b>	<b>8683G</b>	<b>8665G</b>				
 FS-1 24mm	45cm	blue		<b>8629H*</b>	<b>8684G*</b>				

\* New Zealand Indent

^ Australia & New Zealand Indent

# PROLENE

## Monofilament Polypropylene Suture





Needle length and description		Length	Colour	Gauge						
▼ 3/8 CIRCLE REVERSE CUTTING				5/0 1	4/0 1.5	3/0 2	2/0 3	0 3.5	1 4	2 5
	26mm	45cm	blue				<b>8685H</b>			
	30mm	75cm	blue			<b>8675H</b>				
	40mm	75cm	blue			<b>8649H*</b>	<b>8689H</b>			
	75mm	3x50cm	blue							<b>3846T*</b>
▼ 1/2 CIRCLE PRECISION POINT REVERSE CUTTING				5/0 1	4/0 1.5	3/0 2	2/0 3	0 3.5	1 4	2 5
	13mm	45cm	blue	<b>8655G^</b>						
	16mm	45cm	undyed		<b>8603G</b>	<b>DC465^</b>				
▼ STRAIGHT CUTTING				5/0 1	4/0 1.5	3/0 2	2/0 3	0 3.5	1 4	2 5
	60mm	75cm	blue		<b>8621H*</b>	<b>8622H</b>	<b>8623H</b>			
▲ 3/8 CIRCLE CONVENTIONAL CUTTING P NEEDLE				5/0 1	4/0 1.5	3/0 2	2/0 3	0 3.5	1 4	2 5
	26mm	45cm	blue		<b>W8020T^</b>					
▲ 3/8 CIRCLE CONVENTIONAL CUTTING PRIME				6/0 0.7	5/0 1	4/0 1.5	3/0 2	2/0 3	0 3.5	1 4
	13mm	45cm	blue	<b>8617G</b>	<b>8618G*</b>					
	16mm	45cm	blue			<b>W525*</b>				



\* New Zealand Indent











^ Australia & New Zealand Indent






# PROLENE

## Monofilament Polypropylene Suture

Needle length and description		Length	Colour	Gauge						
<b>▲ 3/8 CIRCLE CONVENTIONAL CUTTING PRIME</b>				6/0 0.7	5/0 1	4/0 1.5	3/0 2	2/0 3	0 3.5	1 4
SLIM BLADE NEEDLE	PC-3 	16mm	45cm	blue	W523*	W527*				
	PC-3 	16mm	45cm 45cm	blue blue	W8005T^ 8636G*	W8006T^ 8635G	W8007T^			
	PC-5 	19mm	45cm 45cm	blue blue		8630G^	8631G W8011T^	8632G		
	PC-25 	26mm	45cm	blue				W538*		

<b>▲ 3/8 CIRCLE CONVENTIONAL CUTTING</b>				6/0 0.7	5/0 1	4/0 1.5	3/0 2	2/0 3	0 3.5	1 4
	CFS-3 	16mm	75cm	blue			W8549^			
	SB-35 	36mm	75cm	blue			W8626^			

<b>◎ 3/8 CIRCLE TAPER POINT EVERPOINT™</b>				9/0 0.3	8/0 0.4	7/0 0.5	6/0 0.7	5/0 1	4/0 1.5	3/0 2
	BV130-5 	6.5mm 6.5mm	45cm 60cm	blue blue		EP8730H EP8732H				
	BV175-6 	8mm 8mm	60cm	blue		EP8741H	EP8735H			
	BV175-7 	8mm 8mm	60cm 2 x 60cm	blue blue		EP8767H	EP8766H EPM8766			
	BV175-8 	9.3mm 9.3mm	60cm 75cm	blue blue			EP8747H EP8755H			
	BV-1 	9.3mm 9.3mm	60cm 75cm	blue blue			EP8702H EP8703H	EP8805H EP8709H		
HEMOSEAL	BV-1 	9.3mm 9.3mm	45cm	blue				EPH9702H		
	C-1 	13mm	75cm	blue			EP8889H			
	C-1 	13mm 13mm	45cm 75cm 4x75cm	blue blue blue			EP8800H EP8708H	EP8706H EPM8706		
HEMOSEAL	C-1 	13mm	75cm	blue				EPH8890H		
HEMOSEAL	C-1 	13mm 13mm	60cm 75cm 90cm	blue blue blue				EPH8725H EPH7477H EPH8720H		























<b>◎ 3/8 CIRCLE TAPER POINT</b>				10/0 0.2	9/0 0.3	8/0 0.4	7/0 0.5	6/0 0.7	5/0 1	4/0 1.5
	BV100-4 	5.1mm	15cm	blue	2794G^	2793G^				
	BV100-4 	5.1mm 5.1mm	5cm	blue	2795G^					
	BV175-5 	6.5mm 6.5mm	60cm	blue				D9216^		
	BV130-5 	6.5mm 6.5mm	60cm	blue			8732H*			
	BV130-5 	6.5mm 6.5mm	45cm	blue			8730H			

\* New Zealand Indent

^ Australia & New Zealand Indent

# PROLENE

## Monofilament Polypropylene Suture


















Needle length and description			Length	Colour	Gauge						
⊙ 3/8 CIRCLE TAPER POINT					10/0 0.2	9/0 0.3	8/0 0.4	7/0 0.5	6/0 0.7	5/0 1	4/0 1.5
	BV130-5 	6.5mm 6.5mm	4x45cm 2x60cm 4x60cm	blue			M8730^ M8739^ M8732^				
	BV175-6 	8mm 8mm	60cm	blue			8741H	8735H			
VISIBLACK	BV175-6 	8mm 8mm	60cm	blue			8841H*				
	BV175-6 	8mm 8mm	4x60cm	blue				M8735^			
	BV175-6 	8mm 8mm	75cm	blue			D8292^	D8155^			
	BV175-7 	8mm 8mm	2x45cm	blue				M8765^			
	BV175-7 	8mm 8mm	4x45cm	blue				M8764*			
VISIBLACK	BV175-7 	8mm 8mm	60cm	blue				8966H*			
	BV175-7 	8mm 8mm	60cm	blue			8767H*	8766H*			
	BV175-7 	8mm 8mm	2x60cm	blue			M8767^	M8766*			
PLEDGETS: SOFT 1/8" x 1/8"	BV-1 	9.3mm	4x60cm	blue						D9342^	
PLEDGETS: SOFT 3/16" x 1/8"	BV-1 	9.3mm	4x60cm	blue					D6619^		
	BV-1 	9.3mm 9.3mm	45cm 45cm	blue blue					8806H^ W8712^		
	BV-1 	9.3mm 9.3mm	2x45cm 4x45cm	blue blue				D9511^			
	BV-1 	9.3mm 9.3mm	60cm	blue				8702H	8805H		
	BV-1 	9.3mm 9.3mm	4x60cm	blue				M8702	M8805*		
HEMOSEAL	BV-1 	9.3mm 9.3mm	60cm	blue						9702H	
VISIBLACK	BV-1 	9.3mm 9.3mm	60cm	blue				8304H	8305H^		
VISIBLACK	BV-1 	9.3mm 9.3mm	4x60cm	blue				M8304*			
	BV-1 	9.3mm 9.3mm	75cm	blue				8703H	8709H*		
	BV-1 	9.3mm 9.3mm	4x75cm	blue				M8703^			
	BV175-8 	9.3mm 9.3mm	45cm 60cm	blue blue				8746H* 8747H^			

\* New Zealand Indent

^ Australia & New Zealand Indent

# PROLENE

## Monofilament Polypropylene Suture



Needle length and description			Length	Colour	Gauge						
⊙ 3/8 CIRCLE TAPER POINT					10/0 0.2	9/0 0.3	8/0 0.4	7/0 0.5	6/0 0.7	5/0 1	4/0 1.5
VISIBLACK	BV175-8 	9.3mm 9.3mm	45cm	blue				EH7225H <sup>^</sup>			
	BV175-8 	9.3mm 9.3mm	45cm	blue			8752H <sup>^</sup>				
	BV175-8 	9.3mm 9.3mm	2x45cm	blue			M8752 <sup>^</sup>				
VISIBLACK	BV175-8 	9.3mm 9.3mm	75cm	blue				EH7222H			
	BV175-8 	9.3mm 9.3mm	2x60cm	blue				M8747*			
	BV175-8 	9.3mm 9.3mm	75cm	blue			8754H*	8755H*			
	BV175-8 	9.3mm 9.3mm	2x75cm	blue			M8754 <sup>^</sup>	M8755 <sup>^</sup>			
	BV 	10.5mm 10.5mm	60cm	blue					W8597 <sup>^</sup>		
	C-1 	13mm	75cm	blue					8889H	8890H	
	C-1 	13mm 13mm	45cm	blue				8800H*	8718H*		
	C-1 	13mm 13mm	4x45cm	blue				M8800 <sup>^</sup>	M8718 <sup>^</sup>		
	C-1 	13mm 13mm	60cm 75cm	blue blue				8708H <sup>^</sup>	W8718 <sup>^</sup> 8706H		
	C-1 	13mm 13mm	4x75cm	blue					M8706		
	C-1 	13mm 13mm	90cm	blue					D4946 <sup>^</sup>	8720H	
	C-1 	13mm 13mm	4x90cm	blue						M8720 <sup>^</sup>	
HEMOSEAL	C-1 	13mm 13mm	75cm	blue							9706H <sup>^</sup>
VISIBLACK	C-1 	13mm 13mm	75cm	blue					8307H		



\* New Zealand Indent




<sup>^</sup> Australia & New Zealand Indent







# PROLENE

## Monofilament Polypropylene Suture

Needle length and description			Length	Colour	Gauge						
<b>3/8 CIRCLE TAPER POINT</b>					10/0 0.2	9/0 0.3	8/0 0.4	7/0 0.5	6/0 0.7	5/0 1	4/0 1.5
VISIBLACK	C-1 	13mm 13mm	4x75cm	blue					M8307^		
	C-1 	13mm 13mm	90cm	blue						8321H*	
VISIBLACK	C-1 	13mm 13mm	4x90cm	blue						M8321^	
	BB 	17mm	75cm	blue							8881H
	BB 	17mm 17mm	60cm 90cm	blue blue						8860H^ 8580H^	8581H

<b>1/2 CIRCLE TAPER POINT EVERPOINT™</b>					10/0 0.2	9/0 0.3	8/0 0.4	7/0 0.5	6/0 0.7	5/0 1	4/0 1.5
	RB-2 	13mm 13mm	75cm	blue					EP8711H		
HEMOSEAL	RB-2 	13mm 13mm	75cm 90cm	blue blue						EPH8710H EPH8716H	

<b>1/2 CIRCLE TAPER POINT</b>					10/0 0.2	9/0 0.3	8/0 0.4	7/0 0.5	6/0 0.7	5/0 1	4/0 1.5
	RB-3 	10.5mm	60cm	blue				8206H*			
	RB-2 	13mm 13mm	75cm	blue					8711H		
	RB-2 	13mm 13mm	4x75cm	blue						M8710^	

<b>1/2 CIRCLE TAPER POINT</b>					5/0 1	4/0 1.5	3/0 2	2/0 3	0 3.5	1 4	2 5
	TF 	13mm 13mm	60cm 90cm	blue blue	8205H	8204H D4295^					
	RB-1 	17mm	75cm	blue	8870H	8871H	8872H^				
	RB-1 	17mm 17mm	60cm 90cm 120cm	blue blue blue blue	8555H 8556H W8556^	8557H W8557^ D7377	8558H	8559H^			
	RB-1 	17mm 17mm	4x90cm	blue		M8557^					
HEMOSEAL	RB-1 	17mm 17mm	75cm 90cm	blue blue	HS6856H HS8556H	HS6857H	9556H^				
	SH-2 	20mm 20mm	90cm	blue		W8761					














\* New Zealand Indent

^ Australia & New Zealand Indent



# PROLENE

## Monofilament Polypropylene Suture











Needle length and description		Length	Colour	Gauge							
				5/0 1	4/0 1.5	3/0 2	2/0 3	0 3.5	1 4	2 5	
<b>⊙ 1/2 CIRCLE TAPER POINT</b>											
HEMOSEAL	SH-2 	20mm 20mm	90cm	blue		HS6861H					
	SH-1 	22mm 22mm	75cm	blue			8762H				
	MO-6 	26mm	60cm 75cm	blue blue					8418H	D6731^	
	SH 	26mm	75cm 150cm	blue blue		8831H	8832H	8833H D8336^	8834H		
	SH 	26mm 26mm	90cm	blue			P8522*				
PLEDGETS: FIRM 1/4" x 1/8" x 1/16"											
HEMOSEAL	SH 	26mm 26mm	90cm	blue		HS6855H	HS6822H				
	SH 	26mm 26mm	90cm 90cm 120cm 120cm 150cm	blue blue blue blue blue		8521H D7160^	8522H W8522^ 8534H W8534^	8523H 8533H D8336			
	CT-2 	27mm	75cm 150cm	blue blue				8411H	8412H	D7580^	
	CT-1 	36mm	75cm	blue			8422H*	8423H	8424H	8425H	
	MH 	36mm 36mm	90cm	blue			8842H*	8843H			
	CT 	40mm	75cm	blue					8434H*	8435H	
	CTX 	48mm	75cm	blue						8455H*	
	TP-1 	65mm	150cm 150cm loop	blue blue					D5093		8825G^

\* New Zealand Indent

^ Australia & New Zealand Indent

# PROLENE

## Monofilament Polypropylene Suture






Needle length and description		Length	Colour	Gauge							
				5/0 1	4/0 1.5	3/0 2	2/0 3	0 3.5	1 4	2 5	
<b>1/2 CIRCLE TAPER POINT</b>											
 XLH 70mm 70mm		75cm	blue							8190G <sup>^</sup>	
 XLH 70mm		150cm	blue							8845G <sup>^</sup>	
<b>1/4 CIRCLE TAPERCUT</b>				10/0 0.2	9/0 0.3	8/0 0.4	7/0 0.5	6/0 0.7	5/0 1	4/0 1.5	
 CIF-4 13mm 13mm		20cm	blue	788G <sup>*</sup>							
<b>3/8 CIRCLE TAPERCUT EVERPOINT™</b>				10/0 0.2	9/0 0.3	8/0 0.4	7/0 0.5	6/0 0.7	5/0 1	4/0 1.5	
 CC 9.3mm 9.3mm		60cm	blue				EP8704H	EP8807H			
 CC-1 13mm 13mm		45cm 75cm 4x75cm	blue blue blue					EP8719H EP8707H EPM8707			
<b>3/8 CIRCLE TAPERCUT</b>				10/0 0.2	9/0 0.3	8/0 0.4	7/0 0.5	6/0 0.7	5/0 1	4/0 1.5	
 CC 9.3mm 9.3mm		45cm	blue				8705H <sup>*</sup>				
 CC 9.3mm 9.3mm		4x45cm	blue				M8705 <sup>*</sup>				
 CC 9.3mm 9.3mm		60cm	blue				8704H <sup>*</sup>	8807H <sup>^</sup>	W8803 <sup>^</sup>		
 CC 9.3mm 9.3mm		4x60cm	blue				M8704 <sup>*</sup>	M8807 <sup>^</sup>			
 CC 9.3mm 9.3mm		60cm 75cm	blue blue						W8664 <sup>^</sup>	W8665 <sup>^</sup>	
HEMOSEAL											






\* New Zealand Indent


<sup>^</sup> Australia & New Zealand Indent


# PROLENE


## Monofilament Polypropylene Suture


Needle length and description		Length	Colour	Gauge						
<b>3/8 CIRCLE TAPERCUT</b>				10/0 0.2	9/0 0.3	8/0 0.4	7/0 0.5	6/0 0.7	5/0 1	4/0 1.5
	13mm 13mm	4x75cm	blue					<b>M8707<sup>^</sup></b>		
	13mm 13mm	45cm	blue					<b>8719H<sup>*</sup></b>		
	13mm 13mm	4x45cm	blue					<b>M8719<sup>*</sup></b>		
	13mm 13mm	75cm 90cm	blue blue					<b>8707H<sup>*</sup></b>	<b>8721H<sup>*</sup></b>	
	13mm 13mm	4x90cm	blue						<b>M8721<sup>^</sup></b>	

<b>1/2 CIRCLE TAPERCUT</b>				5/0 1	4/0 1.5	3/0 2	2/0 3	0 3.5	1 4	2 5
	17mm 17mm	90cm	blue		<b>P8935<sup>*</sup></b>					
<b>PLEDGETS: FIRM 1/4" x 1/8" x 1/16"</b>										
	17mm 17mm	90cm	blue	<b>8934H</b>	<b>8935H</b>	<b>8936H<sup>*</sup></b>	<b>8937H<sup>^</sup></b>			
	26mm 26mm	90cm	blue		<b>8975H</b>	<b>8976H</b>	<b>8977H</b>			
	31mm	100cm	blue					<b>W982<sup>*</sup></b>		
	36mm	75cm	blue					<b>8444H</b>		

<b>3/8 CIRCLE TAPERCUT</b>				5/0 1	4/0 1.5	3/0 2	2/0 3	0 3.5	1 4	2 5
	17mm 17mm	90cm	blue			<b>8926H<sup>^</sup></b>				

<b>STRAIGHT TAPER POINT</b>				5/0 1	4/0 1.5	3/0 2	2/0 3	0 3.5	1 4	2 5
	65mm 65mm	75cm	blue					<b>D10061</b>		

<b>MICROPOINT SPATULA STRAIGHT TRANSCAMBER</b>				10/0 0.2	9/0 0.3	8/0 0.4	7/0 0.5	6/0 0.7	5/0 1	4/0 1.5
	16mm 16mm	20cm 23cm	blue blue	<b>1713G</b>	<b>W1713<sup>^</sup></b>					





<b>3/8 CIRCLE MICROPOINT SPATULA</b>				10/0 0.2	9/0 0.3	8/0 0.4	7/0 0.5	6/0 0.7	5/0 1	4/0 1.5
	6.5mm 6.5mm	15cm	blue	<b>1757G<sup>^</sup></b>						

\* New Zealand Indent

<sup>^</sup> Australia & New Zealand Indent

# PROLENE

## Monofilament Polypropylene Suture

Needle length and description		Length	Colour	Gauge						
 <b>1/2 CIRCLE ULTIMA SPATULA</b>				10/0 0.2	9/0 0.3	8/0 0.4	7/0 0.5	6/0 0.7	5/0 1	4/0 1.5
CS160-6 	5.5mm 5.5mm	30cm	blue	<b>7761G*</b>						
 <b>3/8 CIRCLE ULTIMA SPATULA</b>				10/0 0.2	9/0 0.3	8/0 0.4	7/0 0.5	6/0 0.7	5/0 1	4/0 1.5
CS140-6 	6.5mm 6.5mm	30cm	blue	<b>1710G^</b>						

\* New Zealand Indent

^ Australia & New Zealand Indent