

The background is a solid, vibrant blue. Overlaid on this are several thin, white, semi-transparent lines that appear to be made of a light fabric or paper. These lines are draped and folded, creating a sense of movement and depth. They originate from the top left and bottom center, extending towards the right side of the frame. The lines overlap each other, creating varying shades of light blue where they intersect.











**PRONOVA°**

## **PRONOVA\* Suture Poly (hexafluoropropylene-VDF) Suture**

<b>Product Characteristics</b>	<b>PRONOVA* Suture</b>
Physical Material	Polymer of polypropylene
Construction	Monofilament
Absorbable / Non-absorbable	Non-absorbable
Coating (if applicable)	N/A
Colours	Clear and blue
Available Size Range	5/0 through 8/0
Tensile Strength (Wound Support)	Long term wound support
Tensile Strength Retention in Tissue	Indefinite
Absorption Time	Non-absorbable
Frequent Uses	General soft tissue approximation and/or ligation including use in Cardiovascular and Vascular surgery as well as Ophthalmic and Neurological procedures

\*Trademark

# PRONOVA (Hexafluoropropylene-VDF) Monofilament Suture

Needle length and description		Length	Colour	Gauge						
				11/0 0.1	10/0 0.2	9/0 0.3	8/0 0.4	7/0 0.5	6/0 0.7	5/0 1
<b>3/8 CIRCLE TAPER POINT</b>										
 BV130-5	6.5mm	45cm	blue				<b>3730H^</b>			
	6.5mm	60cm	blue				<b>3732H^</b>			
 BV-1 <small>VISIBLACK</small>	9.3mm	60cm	blue					<b>3304H^</b>		
	9.3mm									
 BV-1	9.3mm	60cm	blue					<b>3703H^</b>	<b>3805H*</b>	
	9.3mm	75cm	blue							
 C-1	13mm	75cm	blue						<b>3706H^</b>	
	13mm									
<b>1/2 CIRCLE TAPER POINT</b>										
 RB-1	17mm	60cm	blue	<b>3555H^</b>						
	17mm	90cm	blue	<b>3556H^</b>						
<b>3/8 CIRCLE PRECISION POINT REVERSE CUTTING</b>										
 P-1	11mm	45cm	blue	<b>PN489G</b>						
 P-3	13mm	45cm	blue	<b>PN462G</b>	<b>D9744</b>					
 PS-3	16mm	45cm	blue	<b>PN498G</b>	<b>PN500G</b>	<b>PN501G</b>				
 PS-2	19mm	45cm	blue		<b>PN512G</b>	<b>PN496G</b>				
 PS	26mm	45cm	blue				<b>PN584G</b>			

\* New Zealand Indent

^ Australia & New Zealand Indent

