

The background features a solid light purple color with several translucent, white, flowing lines that create a sense of movement and depth, resembling liquid or smoke. These lines originate from the top left and bottom center, sweeping across the frame.

**Coated
VICRYL® PLUS**

Coated VICRYL* Plus Antibacterial (polyglactin 910) Suture

Product Characteristics	Coated VICRYL* Plus Suture
Physical Material	Polyglactin 910 (copolymer of glycolide and lactide)
Construction	Braided (Monofilament in sizes 9/0 and 10/0)
Absorbable / Non-absorbable	Absorbable
Coating (if applicable)	50% Polyglactin 370 and 50% Calcium Stearate
Colours	Undyed beige and Dyed violet
Available Size Range	2 through 5/0
Tensile Strength (Wound Support)	Medium Term (4 weeks)
Tensile Strength Retention in Tissue	75% of original strength remains at 2 weeks 50% of original strength remains at 3 weeks 25% of original strength remains at 4 weeks
Absorption Time	Absorbed in 56-70 days, 63 days average
Frequent Uses	General soft tissue approximation and/or ligation, including use in Ophthalmic procedures

VICRYL™ Plus Suture with IRGACARE® MP (triclosan) are effective against the pathogens most commonly associated with Surgical site infections:^{1,2,3}

- *Staphylococcus aureus*
- *Staphylococcus epidermidis*
- Methicillin-resistant *Staphylococcus aureus* (MRSA)
- Methicillin-resistant *Staphylococcus aureusepidermidis* (MRSE)












*Trademark

¹ Wang ZX, Jiang CP, Cao Y, Ding DT. Systematic Review and meta-analysis of triclosan-coated sutures for the prevention of surgical site infection. *British Journal of Surgery* 2013;100:465-473















² Edmiston CE, Daoud FC, Leaper D. Is there an evidence-based argument for embracing an antimicrobial (triclosan) coated suture technology to reduce the risk for surgical site infections?: A meta-analysis. *Surgery* 2103:154;89-100

³ Coated VICRYL Plus IFU (instructions for use)












Coated VICRYL® Plus Antibacterial Suture (Polyglactin 910) Braided Suture

Needle length and description		Length	Colour	Gauge						
○ 1/2 CIRCLE ETHIGUARD BLUNT POINT				3/0 2	2/0 3	0 3.5	1 4	2 5	3 6	4 6
SHB 	26mm	8x45cm control release	undyed							
CTB-2 	26mm	70cm	violet							
MHB-1 	31mm	75cm	violet							
CTB-1 	36mm	8x45cm control release	undyed			VCPB840D^				
CTB-1 	36mm	70cm 70cm 75cm 90cm 90cm	violet undyed violet undyed			VCPB340H^	VCPB341H^			
					VCPB945H*	VCPB946H				
CTB 	40mm	75cm 90cm	violet violet							
BT-3 	45mm	90cm	violet							
CTX-B 	48mm	8x45cm control release	undyed			VCPB724D^	VCPB725D^			
CTX-B 	48mm	90cm 90cm	violet undyed			VCPB978H*	VCPB977H			
BP-1 	65mm	4x90cm	violet							
BP-1 	65mm	100cm	violet							


Coated VICRYL® Plus Antibacterial Suture (Polyglactin 910) Braided Suture




Needle length and description		Length	Colour	Gauge						
▼ COMPUND CURVE PRECISION POINT REVERSE CUTTING				5/0 1	4/0 1.5	3/0 2	2/0 3	0 3.5	1 4	2 5
 PS-4C	16mm	45cm	undyed		VCP656G					
▼ 3/8 CIRCLE PRECISION POINT REVERSE CUTTING				7/0 0.5	6/0 0.7	5/0 1	4/0 1.5	3/0 2	2/0 3	0 3.5
 P-1	10.5mm	45cm 45cm	undyed undyed							
 P-3	13mm	45cm	undyed			VCP493G	VCP494G			
 PS-3	16mm	45cm	undyed			VCP500G				
 PS-2	19mm	45cm 45cm	undyed undyed				VCP496H	VCP497H		
 PS-1	24mm	45cm	undyed					VCP683H		
▼ 3/8 CIRCLE REVERSE CUTTING				7/0 0.5	6/0 0.7	5/0 1	4/0 1.5	3/0 2	2/0 3	0 3.5
 FS-2	19mm	45cm 70cm	violet undyed			VCP421H^	VCP422H	VCP393H^ VCP423H		
 FS-1	24mm	70cm 70cm	violet undyed				VCP451H^ VCP441H*	VCP442H	VCP443H	
 FSL	30mm	70cm 70cm	violet undyed					VCP588H*	VCP589H	
▼ 1/2 CIRCLE PRECISION POINT REVERSE CUTTING				7/0 0.5	6/0 0.7	5/0 1	4/0 1.5	3/0 2	2/0 3	0 3.5
 P-2	8mm	45cm	undyed			VCP503G^	VCP504G^			
▼ 1/2 CIRCLE REVERSE CUTTING				3/0 2	2/0 3	0 3.5	1 4	2 5	3 6	4 6
 OS-4	22mm	3x45cm control release	violet							
 OS-4	22mm	70cm 75cm	undyed violet			VCP694H^	VCP695H			
 X-1	22mm	70cm 70cm	violet undyed	VCP460H^ VCP458H^						
 CP-2	26mm	8x45cm 70cm	violet undyed		VCP762D VCP869H	VCP727D VCP870H^	VCP871H			











Coated VICRYL® Plus Antibacterial Suture (Polyglactin 910) Braided Suture

Needle length and description		Length	Colour	Gauge							
▼ 1/2 CIRCLE REVERSE CUTTING				3/0 2	2/0 3	0 3.5	1 4	2 5	3 6	4 6	
	CP-1	36mm	70cm 70cm 90cm	violet undyed violet		VCP466H^ VCP472H^	VCP467H^ VCP267H^	VCP468H* VCP268H^ VCP474H			
	OS-6	36mm	70cm	undyed				VCP535H			
	CP	40mm	70cm 90cm	violet violet		VCP478H*		VCP480H			
	OS-8	40mm	70cm	undyed							
	CPX	48mm	70cm	undyed				VCP569H			
▼ STRAIGHT CUTTING				4/0 1.5	3/0 2	2/0 3	0 3.5	1 4	2 5	3 6	
	KS	60mm	70cm 70cm 75cm	violet undyed undyed	VCP662H^	VCP663H					
▲ 3/8 CIRCLE CONVENTIONAL CUTTING P NEEDLE				4/0 1.5	3/0 2	2/0 3	0 3.5	1 4	2 5	3 6	
	PC-26	26mm	45cm 75cm	undyed undyed							
▲ 3/8 CIRCLE CONVENTIONAL CUTTING PRIME				6/0 0.7	5/0 1	4/0 1.5	3/0 2	2/0 3	0 3.5	1 4	
	PC-1	13mm	45cm	undyed		VCP843G*	VCP835G*				
	PC-3	16mm	45cm 45cm	undyed undyed		VCP844G*					
	PC-5	19mm	45cm 45cm 75cm	undyed undyed undyed			VCP823G	VCP824G			
▲ 1/2 CIRCLE CONVENTIONAL CUTTING				6/0 0.7	5/0 1	4/0 1.5	3/0 2	2/0 3	0 3.5	1 4	
	CX-1	22mm	45cm	violet							













Coated VICRYL® Plus Antibacterial Suture (Polyglactin 910) Braided Suture

Needle length and description	Length	Colour	Gauge						
▲ 1/2 CIRCLE CONVENTIONAL CUTTING STERNUM			0 3.5	1 4	2 5	3 6	4 6	5 7	7 9
CCS-1 	60mm	70cm	undyed						






◎ 3/8 CIRCLE TAPER POINT				9/0 0.3	8/0 0.4	7/0 0.5	6/0 0.7	5/0 1	4/0 1.5	3/0 2
BV100-3 	3.9mm	15cm	violet							
BV130-5 	6.5mm	15cm	violet							
BV-1 	9.3mm	30cm	violet							

◎ 1/2 CIRCLE TAPER POINT				6/0 0.7	5/0 1	4/0 1.5	3/0 2	2/0 3	0 3.5	1 4
RB-2 	13mm	45cm	undyed							
TF 	13mm	70cm	undyed		VCP433H					
RB-1 	17mm	70cm 70cm	violet undyed		VCP303H	VCP304H VCP214H	VCP305H	VCP306H		
CT-3 	22mm	70cm	violet				VCP327H^		VCP329H^	
SH-1 	22mm	70cm	violet			VCP310H*	VCP311H			
SH 	26mm	8x45cm control release 45cm	violet undyed				VCP774D VCP864D^	VCP775D^		
SH 	26mm	8x70cm control release	violet							
SH 	26mm	70cm 70cm 75cm	violet undyed violet			VCP315H	VCP316H VCP416H	VCP317H VCP417H	VCP418H^	
SH 	26mm 26mm	90cm	violet				VCP527H^			
CT-2 	27mm	8x45cm control release	violet							
























Coated VICRYL® Plus Antibacterial Suture (Polyglactin 910) Braided Suture

Needle length and description		Length	Colour	Gauge						
⊙ 1/2 CIRCLE TAPER POINT				5/0 1	4/0 1.5	3/0 2	2/0 3	0 3.5	1 4	2 5
CT-2 	27mm	70cm 70cm	violet undyed			VCP332H^ VCP232H^	VCP333H VCP269H	VCP334H VCP270H*	VCP335H	
MH-1 	31mm	75cm	violet							
MO-5 	31mm	75cm	violet							
CT-1 	36mm	8x45cm control release	undyed					VCP840D		
CT-1 	36mm	70cm 70cm 90cm 90cm	violet undyed violet undyed			VCP338H^	VCP339H VCP259H VCP345H VCP945H	VCP340H VCP260H VCP346H VCP946H*	VCP341H VCP261H* VCP347H	
MH 	36mm	70cm 75cm	violet violet			VCP322H^	VCP323H			
MO-4 	36mm	70cm	violet					VCP436H*	VCP437H*	
CT 	40mm	8x45cm control release	violet				VCP751D*	VCP752D^	VCP753D^	
CT 	40mm	70cm 70cm 90cm 90cm 90cm	violet undyed violet undyed violet				VCP351H VCP275H* VCP357H*	VCP352H VCP358H VCP958H^	VCP353H VCP359H VCP959H^	
CTX 	48mm	70cm 90cm	violet violet				VCP363H* VCP369H^	VCP364H VCP370H^	VCP365H VCP371H	
TP-1 	65mm	70cm 135cm	violet undyed							VCP880T*
XLH 	70mm	70cm	violet				VCP581G*		VCP583G^	

Coated VICRYL® Plus Antibacterial Suture (Polyglactin 910) Braided Suture

⊙ 5/8 CIRCLE TAPER POINT				2/0 3	0 3.5	1 4	2 5	3 6	4 6	5 7
	27mm	70cm	violet	VCP602H	VCP603H					
	36mm	70cm	violet	VCP375H^						
	40mm	70cm 75cm	violet violet							
▼ 1/4 CIRCLE TAPERCUT				8/0 0.4	7/0 0.5	6/0 0.7	5/0 1	4/0 1.5	3/0 2	2/0 3
	5mm	25cm	violet							
▼ 1/2 CIRCLE TAPERCUT				3/0 2	2/0 3	0 3.5	1 4	2 5	3 6	4 6
	36mm	75cm 90cm 90cm	violet violet undyed	VCP516H*	VCP517H*	VCP518H	VCP519H			

Coated VICRYL® Plus Antibacterial Suture (Polyglactin 910) Braided Suture

Needle length and description	Length	Colour	Gauge						
 1/2 CIRCLE MICROPOINT SPATULA			9/0 0.3	8/0 0.4	7/0 0.5	6/0 0.7	5/0 1	4/0 1.5	3/0 2
TG160-8 	5.5mm 5.5mm	45cm	violet						
 1/4 CIRCLE SABRELOC CONVENTIONAL SPATULA			9/0 0.3	8/0 0.4	7/0 0.5	6/0 0.7	5/0 1	4/0 1.5	3/0 2
S-29 	8mm 8mm	30cm 45cm	violet violet						
 1/2 CIRCLE SABRELOC CONVENTIONAL SPATULA			9/0 0.3	8/0 0.4	7/0 0.5	6/0 0.7	5/0 1	4/0 1.5	3/0 2
S-28 	8mm 8mm	45cm	violet						
 1/4 CIRCLE SABRELOC SPATULA			9/0 0.3	8/0 0.4	7/0 0.5	6/0 0.7	5/0 1	4/0 1.5	3/0 2
S-14 	8mm 8mm	20cm 45cm 45cm	violet violet undyed				VCP671G*		
 1/2 CIRCLE SABRELOC SPATULA			9/0 0.3	8/0 0.4	7/0 0.5	6/0 0.7	5/0 1	4/0 1.5	3/0 2
S-22 	8mm 8mm	45cm 45cm	violet undyed						
 1/4 CIRCLE SABRELOC CENTRE POINT SPATULA			9/0 0.3	8/0 0.4	7/0 0.5	6/0 0.7	5/0 1	4/0 1.5	3/0 2
RD-1 	8mm 8mm	30cm	violet						
 1/8 CIRCLE ULTIMA SPATULA			10/0 0.2	9/0 0.3	8/0 0.4	7/0 0.5	6/0 0.7	5/0 1	4/0 1.5
CS44-6 	5.5mm	10cm	violet						
 1/4 CIRCLE ULTIMA SPATULA			10/0 0.2	9/0 0.3	8/0 0.4	7/0 0.5	6/0 0.7	5/0 1	4/0 1.5
CS90-6 	7mm	10cm	violet						
 3/8 CIRCLE ULTIMA SPATULA			10/0 0.2	9/0 0.3	8/0 0.4	7/0 0.5	6/0 0.7	5/0 1	4/0 1.5
CS140-6 	6.5mm	10cm	violet						
 3/8 CIRCLE ULTIMA SPATULA			10/0 0.2	9/0 0.3	8/0 0.4	7/0 0.5	6/0 0.7	5/0 1	4/0 1.5
CS140-6 	6.5mm	10cm	violet						
CS140-8 	6.75mm 6.75mm	20cm	violet						
 1/2 CIRCLE ULTIMA SPATULA			10/0 0.2	9/0 0.3	8/0 0.4	7/0 0.5	6/0 0.7	5/0 1	4/0 1.5
CS160-6 	5.5mm 5.5mm	30cm	violet						

Coated VICRYL® Plus Antibacterial Suture (Polyglactin 910) Braided Suture

Needle length and description	Length	Colour	Gauge						
NON-NEEDED			5/0 1	4/0 1.5	3/0 2	2/0 3	0 3.5	1 4	2 5
LIGAPAK*	135cm 135cm 360cm	violet undyed undyed			VCP205G VCP885G	VCP286G	VCP207G	VCP208G	
NON-NEEDED	6x45cm	violet					VCP106G		
NON-NEEDED	12x45cm	violet		VCP103G	VCP104G	VCP105G			
NON-NEEDED	135cm 135cm 150cm 150cm	violet undyed violet undyed				VCP615H VCP607H*	VCP616H	VCP617H	VCP618H