








Coated
VICRYL *RAPIDE*[®]



Coated VICRYL *RAPIDE** (polyglactin 910) Suture


Product Characteristics	Coated VICRYL <i>RAPIDE</i>* Suture
Physical Material	Polyglactin 910 (copolymer of glycolide and lactide)
Construction	Braided
Absorbable / Non-absorbable	Absorbable
Coating (if applicable)	50% Polyglactin 370 and 50% Calcium Stearate
Colours	Undyed – beige and Dyed – violet
Available Size Range	1 through 7/0
Tensile Strength (Wound Support)	Short Term (7-10 days)
Tensile Strength Retention in Tissue	50% of original strength remains at 5 days Suture strength is lost between 10-14 days
Absorption Time	Completely absorbed in 42 days
Frequent Uses	Superficial soft tissue approximation of the skin and mucosa where short term wound support is required


*Trademark




Coated VICRYL RAPIDE® (Polyglactin 910) Braided Suture


Needle length and description		Length	Colour	Gauge						
▼ 3/8 CIRCLE PRECISION POINT REVERSE CUTTING				6/0 0.7	5/0 1	4/0 1.5	3/0 2	2/0 3	0 3.5	1 4
 P-1	11mm	45cm 45cm 45cm	violet undyed undyed	V32H W9913	MPVR9915H* W9915					
 P-3	13mm	45cm	undyed		MPVR4930H	MPVR4940H*				
 PS-2	19mm	45cm 75cm 75cm	undyed undyed undyed			MPVR4960H MPVR4260H W9922	MPVR4970H W9923			
 PS-1	24mm	70cm	undyed				MPVR9360H*			
 PS	26mm	75cm	undyed				W9932			

▼ 3/8 CIRCLE REVERSE CUTTING				6/0 0.7	5/0 1	4/0 1.5	3/0 2	2/0 3	0 3.5	1 4
 FS-3	16mm	75cm	undyed		VR2296^					
 FS-1	24mm	45cm	undyed				V4421H^			

▲ 3/8 CIRCLE CONVENTIONAL CUTTING PRIME				6/0 0.7	5/0 1	4/0 1.5	3/0 2	2/0 3	0 3.5	1 4
 PC-3	16mm	45cm	undyed			W9918				

▲ 1/2 CIRCLE CONVENTIONAL CUTTING P NEEDLE				6/0 0.7	5/0 1	4/0 1.5	3/0 2	2/0 3	0 3.5	1 4
 PC-22	22mm	45cm	undyed				W9927			






◎ 1/2 CIRCLE TAPER POINT				6/0 0.7	5/0 1	4/0 1.5	3/0 2	2/0 3	0 3.5	1 4
 RB-1	17mm	70cm	undyed			W9970				
 SH	26mm	70cm	undyed				V4160H	V4170H*		
 CT-1	36mm	90cm	undyed					V9450H		

◐ 1/2 CIRCLE TAPERCUT				6/0 0.7	5/0 1	4/0 1.5	3/0 2	2/0 3	0 3.5	1 4
 V-34	36mm	90cm	undyed					W9962		

* New Zealand Indent

^ Australia & New Zealand Indent

Coated VICRYL *RAPIDE*[®] (Polyglactin 910) Braided Suture

Needle length and description		Length	Colour	Gauge						
 3/8 CIRCLE MICROPPOINT SPATULA				7/0 0.5	6/0 0.7	5/0 1	4/0 1.5	3/0 2	2/0 3	0 3.5
TG140-8 	6.5mm	30cm	violet	W9914*						
 1/2 CIRCLE TAPER POINT, 1/2 CIRCLE REVERSE CUTTING				7/0 0.5	6/0 0.7	5/0 1	4/0 1.5	3/0 2	2/0 3	0 3.5
CT-1 / CP-1 	36mm 36mm	135cm	undyed						VC532H*	
STRAIGHT CUTTING				7/0 0.5	6/0 0.7	5/0 1	4/0 1.5	3/0 2	2/0 3	0 3.5
SC-1 	13mm	45cm	undyed				W9909			

* New Zealand Indent

^ Australia & New Zealand Indent