

The background features a solid light purple color with several translucent, white, flowing lines that create a sense of movement and depth. These lines originate from the top left and bottom center, sweeping across the frame in various directions.













**Coated  
VICRYL®**

## Coated VICRYL\* (polyglactin 910) Suture

Product Characteristics	Coated VICRYL* Suture
Physical Material	Polyglactin 910 (copolymer of glycolide and lactide)
Construction	Braided (Monofilament in sizes 9/0 and 10/0)
Absorbable / Non-absorbable	Absorbable
Coating (if applicable)	50% Polyglactin 370 and 50% Calcium Stearate
Colours	Undyed – beige and Dyed – violet
Available Size Range	3 through 8/0
Tensile Strength (Wound Support)	Medium Term (4 weeks)
Tensile Strength Retention in Tissue	6/0 and larger: 75% of original strength remains at 2 weeks 50% of original strength remains at 3 weeks 25% of original strength remains at 4 weeks  7/0 and smaller: 75% of original strength remains at 2 weeks 40% of original strength remains at 3 weeks
Absorption Time	Absorbed in 56-70 days, 63 days average
Frequent Uses	General soft tissue approximation and/or ligation, including use in Ophthalmic procedures

\*Trademark


# Coated VICRYL® (Polyglactin 910) Braided Suture







Needle length and description		Length	Colour	Gauge						
○ 1/2 CIRCLE ETHIGUARD BLUNT POINT				3/0 2	2/0 3	0 3.5	1 4	2 5	3 6	4 6
SHB 	26mm	8x45cm control release	undyed	JB864^						
CTB-2 	26mm	70cm	violet				D9754*			
MHD-1 	31mm	75cm	violet		W9984*					
CTB-1 	36mm	8x45cm control release	undyed			JB840^	JB841^			
CTB-1 	36mm	70cm 70cm 90cm	violet undyed undyed	JB944^	JB259 JB945	JB340* JB946^	JB341 JB947			
CTD-1 	36mm	75cm 90cm	violet violet		W9993^	W9994^	W9995			
CTD 	40mm	75cm 90cm	violet violet			W9988*	W9997			
TD-3 	45mm	90cm	violet				W9998	W9999*		
CTX-B 	48mm	8x45cm control release	undyed			JB724^	JB725^			
CTX-B 	48mm	90cm 90cm	violet undyed			JB978*	JB977			
BP-1 	65mm	4x90cm	violet					D8659^		
○ 3/8 CIRCLE ETHIGUARD BLUNT POINT				3/0 2	2/0 3	0 3.5	1 4	2 5	3 6	4 6
BP-5 	65mm	100cm	violet				W9391^			




\* New Zealand Indent


^ Australia & New Zealand Indent





# Coated VICRYL® (Polyglactin 910) Braided Suture

Needle length and description		Length	Colour	Gauge						
<b>▼ COMPUND CURVE PRECISION POINT REVERSE CUTTING</b>				5/0 1	4/0 1.5	3/0 2	2/0 3	0 3.5	1 4	2 5
 PS-4C	16mm	45cm	undyed	<b>W9436</b>	<b>J656G</b>					

<b>▼ 3/8 CIRCLE PRECISION POINT REVERSE CUTTING</b>				7/0 0.5	6/0 0.7	5/0 1	4/0 1.5	3/0 2	2/0 3	0 3.5
 P-1	10.5mm	45cm 45cm	undyed undyed	<b>J488G*</b>	<b>J489G</b> <b>W9500T^</b>	<b>W9501T^</b>				
 P-1	11mm 11mm	45cm	violet		<b>W9537</b>					
 P-3	13mm	45cm	undyed		<b>J492G</b>	<b>J493G</b>	<b>J494G</b>			
 PS-3	16mm	45cm	undyed			<b>J500G</b>				
 PS-2	19mm	45cm 45cm	undyed undyed				<b>J496H</b> <b>W9515T^</b>	<b>J497H</b>		
 PS-1	24mm	45cm	undyed					<b>J683H</b>		

<b>▼ 3/8 CIRCLE REVERSE CUTTING</b>				7/0 0.5	6/0 0.7	5/0 1	4/0 1.5	3/0 2	2/0 3	0 3.5
 FS-2	19mm	45cm 70cm	violet undyed			<b>J421H*</b>	<b>J422H</b>	<b>J393H</b> <b>J423H</b>		
 FS-1	24mm	70cm 70cm	violet undyed				<b>J451H^</b> <b>J441H</b>	<b>J442H</b>	<b>J443H</b>	
 FSL	30mm	70cm	undyed					<b>J588H</b>	<b>J589H</b>	












<b>▼ 1/2 CIRCLE PRECISION POINT REVERSE CUTTING</b>				7/0 0.5	6/0 0.7	5/0 1	4/0 1.5	3/0 2	2/0 3	0 3.5
 P-2	8mm	45cm	undyed			<b>J503G^</b>	<b>J504G*</b>			

<b>▼ 1/2 CIRCLE REVERSE CUTTING</b>				3/0 2	2/0 3	0 3.5	1 4	2 5	3 6	4 6
 OS-4	22mm	45cm control release	violet				<b>J708T^</b>			
 OS-4	22mm	70cm 75cm	undyed violet		<b>J693H^</b>	<b>J694H^</b>	<b>J695H</b>	<b>W9287^</b>		
 X-1	22mm	70cm 70cm	violet undyed	<b>J460H^</b> <b>J458H^</b>						
 CP-2	26mm	70cm	undyed		<b>J869H</b>	<b>J870H</b>	<b>J871H</b>			

\* New Zealand Indent

^ Australia & New Zealand Indent















# Coated VICRYL® (Polyglactin 910) Braided Suture

Needle length and description		Length	Colour	Gauge						
<b>▼ 1/2 CIRCLE REVERSE CUTTING</b>				3/0 2	2/0 3	0 3.5	1 4	2 5	3 6	4 6
 CP-1	36mm	70cm 70cm 90cm	violet undyed violet	J471H^	J466H J472H*	J467H* J267H*	J468H* J268H^ J474H^			
 OS-6	36mm	70cm	undyed			J534H^	J535H			
 CP	40mm	70cm 90cm	violet violet		J484H^	J485H^				
 OS-8	40mm	70cm	undyed			J699H^				
 CPX	48mm	70cm	undyed				J569H			
<b>▼ STRAIGHT CUTTING</b>				4/0 1.5	3/0 2	2/0 3	0 3.5	1 4	2 5	3 6
 KS	60mm	70cm 70cm 75cm	violet undyed undyed	J662H*	J663H W9717^	J664H				
<b>▲ 3/8 CIRCLE CONVENTIONAL CUTTING P NEEDLE</b>				4/0 1.5	3/0 2	2/0 3	0 3.5	1 4	2 5	3 6
 PC-26	26mm	45cm 75cm	undyed undyed	W9520T^ W9580T^	W9521T^ W9581T^	W9522T W9582T^				
<b>▲ 3/8 CIRCLE CONVENTIONAL CUTTING PRIME</b>				6/0 0.7	5/0 1	4/0 1.5	3/0 2	2/0 3	0 3.5	1 4
 PC-1	13mm	45cm	undyed	J833G	J834G	J835G				
 PC-3	16mm	45cm 45cm	undyed undyed		J844G* W9505T^	W9506T				
 PC-5	19mm	45cm 45cm 75cm	undyed undyed undyed			J823G W9510T^ W9570T^	J824G* W9511T^ W9571T^			
<b>▲ 1/2 CIRCLE CONVENTIONAL CUTTING</b>				6/0 0.7	5/0 1	4/0 1.5	3/0 2	2/0 3	0 3.5	1 4
 CX-1	22mm	45cm	violet				W9472^			

\* New Zealand Indent

^ Australia & New Zealand Indent














# Coated VICRYL® (Polyglactin 910) Braided Suture

Needle length and description	Length	Colour	Gauge							
<b>▲ 1/2 CIRCLE CONVENTIONAL CUTTING STERNUM</b>			0 3.5	1 4	2 5	3 6	4 6	5 7	7 9	
 <p>CCS-1</p>	60mm	70cm	undyed				<b>J808T*</b>			
<b>◎ 3/8 CIRCLE TAPER POINT</b>			9/0 0.3	8/0 0.4	7/0 0.5	6/0 0.7	5/0 1	4/0 1.5	3/0 2	
 <p>BV100-3</p>	3.9mm	15cm	violet	<b>V402G</b>						
 <p>BV130-5</p>	6.5mm	15cm	violet		<b>J401G</b>					
 <p>BV-1</p>	9.3mm	30cm	violet		<b>W9577^</b>					
<b>◎ 1/2 CIRCLE TAPER POINT</b>			6/0 0.7	5/0 1	4/0 1.5	3/0 2	2/0 3	0 3.5	1 4	
 <p>RB-2</p>	13mm	45cm	undyed	<b>W9981</b>	<b>W9982*</b>					
 <p>TF</p>	13mm	70cm	undyed		<b>J433H</b>					
 <p>RB-1</p>	17mm	70cm 70cm	violet undyed		<b>J303H</b>	<b>J304H J214H</b>	<b>J305H</b>	<b>J306H</b>		
 <p>CT-3</p>	22mm	70cm	violet				<b>J327H^</b>		<b>J329H^</b>	
 <p>SH-1</p>	22mm	70cm	violet			<b>J310H</b>	<b>J311H</b>			
 <p>SH</p>	26mm	8x45cm control release	violet undyed				<b>J774D J864D*</b>	<b>J775D^</b>		
 <p>SH</p>	26mm	70cm control release	violet				<b>J784G^</b>			
 <p>SH</p>	26mm	70cm 70cm 75cm	violet undyed violet			<b>J315H</b>	<b>J316H J416H W9120^</b>	<b>J317H J417H W9121^</b>	<b>J318H</b>	
 <p>SH</p>	26mm 26mm	90cm	violet				<b>J527H^</b>			
 <p>CT-2</p>	27mm	8x45cm control release	violet					<b>J726D^</b>	<b>J727D^</b>	

\* New Zealand Indent

^ Australia & New Zealand Indent








# Coated VICRYL® (Polyglactin 910) Braided Suture

Needle length and description	Length	Colour	Gauge							
◎ 1/2 CIRCLE TAPER POINT			4/0 1.5	3/0 2	2/0 3	0 3.5	1 4	2 5	3 6	
CT-2 	27mm	70cm 70cm	violet undyed		J332H J232H^	J333H J269H	J334H J270H^	J335H		
MH-1 	31mm	75cm	violet					W9213^		
MO-5 	31mm	75cm	violet					W9216*		
CT-1 	36mm	8x45cm control release	undyed				J840D^			
CT-1 	36mm	70cm 70cm 90cm 90cm	violet undyed violet undyed		J338H*	J339H J259H J345H J945H	J340H J260H J346H J946H*	J341H J261H* J347H		
MH 	36mm	70cm 75cm	violet violet		J322H*	J323H	J324H W9141^			
MO-4 	36mm	70cm 75cm	violet undyed				J436H*	J437H*		D8768^
CT 	40mm	8x45cm control release	violet			J751D^	J752D*	J753D^		
CT 	40mm	70cm 70cm 90cm 90cm 90cm	violet undyed violet undyed violet			J351H J275H* J357H*	J352H J358H J958H^	J353H J359H J959H^ W9431^		
MO-45 	45mm	75cm	violet						W9246	
CTX 	48mm	70cm 90cm	violet violet			J363H J369H	J364H J370H^	J365H J371H		
TP-1 	65mm	4x70cm 70cm 135cm 135cm	violet violet undyed violet					J650G^ D5792	J649G J880T*	
XLH 	70mm	70cm	violet			J581G*		J583G*		

\* New Zealand Indent

^ Australia & New Zealand Indent

# Coated VICRYL® (Polyglactin 910) Braided Suture











Needle length and description		Length	Colour	Gauge							
<b>1/2 CIRCLE TAPER POINT, 1/2 CIRCLE REVERSE CUTTING</b>				6/0 0.7	5/0 1	4/0 1.5	3/0 2	2/0 3	0 3.5	1 4	
CT-1 / CP-1 		36mm 36mm	135cm	violet					<b>J532H*</b>		
<b>5/8 CIRCLE TAPER POINT</b>				3/0 2	2/0 3	0.3 5	1 4	2 5	3 6	4 6	
UR-6 		27mm	70cm	violet	<b>D7973^</b>	<b>J602H</b>	<b>J603H</b>				
UR-6 		27mm 27mm	90cm	undyed		<b>D8923</b>					
UR-5 		36mm	70cm	violet		<b>J375H^</b>					
UR-4 		40mm	70cm 75cm	violet violet		<b>J381H^</b> <b>W9165^</b>					
<b>1/4 CIRCLE TAPERCUT</b>				8/0 0.4	7/0 0.5	6/0 0.7	5/0 1	4/0 1.5	3/0 2	2/0 3	
V130-4 		5mm	25cm	violet	<b>J408G*</b>						
<b>1/2 CIRCLE TAPERCUT</b>				3/0 2	2/0 3	0 3.5	1 4	2 5	3 6	4 6	
V-34 		36mm	75cm 90cm 90cm	violet violet undyed	<b>J516H^</b> <b>J916H^</b>	<b>J517H</b>	<b>J518H</b>	<b>W9365^</b> <b>J519H</b>	<b>W9357^</b>		

\* New Zealand Indent

^ Australia & New Zealand Indent
























# Coated VICRYL® (Polyglactin 910) Braided Suture

Needle length and description		Length	Colour	Gauge						
<b>1/2 CIRCLE TAPERCUT</b>				3/0 2	2/0 3	0 3.5	1 4	2 5	3 6	4 6
	V-39 45mm	75cm	violet			W9376*	W9377	W9378		
<b>5/8 CIRCLE TAPERCUT</b>				6/0 0.7	5/0 1	4/0 1.5	3/0 2	2/0 3	0 3.5	1 4
	UV-10 10mm	30cm	violet	W9172*						
<b>J NEEDLE TAPERCUT</b>				6/0 0.7	5/0 1	4/0 1.5	3/0 2	2/0 3	0 3.5	1 4
	V-38 31mm	75cm	violet							W9335*
	J 22mm	75cm	violet				W9223	W9220	W9221	
<b>3/8 CIRCLE MICROPOINT REVERSE CUTTING</b>				6/0 0.7	5/0 1	4/0 1.5	3/0 2	2/0 3	0 3.5	1 4
	G-6 8mm 8mm	45cm	violet	J551G*						
<b>1/4 CIRCLE MICROPOINT SPATULA</b>				6/0 0.7	5/0 1	4/0 1.5	3/0 2	2/0 3	0 3.5	1 4
	TG100-8 6.5mm 6.5mm	45cm	violet	J544G*						
<b>3/8 CIRCLE MICROPOINT SPATULA</b>				9/0 0.3	8/0 0.4	7/0 0.5	6/0 0.7	5/0 1	4/0 1.5	3/0 2
	TG140-6 6.5mm 6.5mm	30cm	violet		J974G*					
	TG140-8 6.5mm 6.5mm	30cm	violet			W9561*				
	TG140-8 6.5mm 6.5mm	30cm 45cm	violet violet	V549G*	J548G	J566G J546G				
<b>1/2 CIRCLE MICROPOINT SPATULA</b>				9/0 0.3	8/0 0.4	7/0 0.5	6/0 0.7	5/0 1	4/0 1.5	3/0 2
	TG160-6 5.5mm 5.5mm	30cm	violet		J975G*					

\* New Zealand Indent

^ Australia & New Zealand Indent

# Coated VICRYL® (Polyglactin 910) Braided Suture

Needle length and description		Length	Colour	Gauge						
 <b>1/2 CIRCLE MICROPOINT SPATULA</b>				9/0 0.3	8/0 0.4	7/0 0.5	6/0 0.7	5/0 1	4/0 1.5	3/0 2
TG160-8 	5.5mm 5.5mm	45cm	violet			<b>J576G</b>				
 <b>1/4 CIRCLE SABRELOC CONVENTIONAL SPATULA</b>				9/0 0.3	8/0 0.4	7/0 0.5	6/0 0.7	5/0 1	4/0 1.5	3/0 2
S-29 	8mm 8mm	30cm 45cm	violet violet				<b>J556G</b> <b>J555G*</b>			
 <b>1/2 CIRCLE SABRELOC CONVENTIONAL SPATULA</b>				9/0 0.3	8/0 0.4	7/0 0.5	6/0 0.7	5/0 1	4/0 1.5	3/0 2
S-28 	8mm 8mm	45cm	violet				<b>J562G</b>			
 <b>1/4 CIRCLE SABRELOC SPATULA</b>				9/0 0.3	8/0 0.4	7/0 0.5	6/0 0.7	5/0 1	4/0 1.5	3/0 2
S-14 	8mm 8mm	20cm 45cm 45cm	violet violet undyed				<b>W9562^</b> <b>J570G</b> <b>J670G</b>	<b>J571G</b> <b>J671G*</b>		
 <b>1/2 CIRCLE SABRELOC SPATULA</b>				9/0 0.3	8/0 0.4	7/0 0.5	6/0 0.7	5/0 1	4/0 1.5	3/0 2
S-22 	8mm 8mm	45cm 45cm	violet undyed				<b>W9566</b> <b>W9756</b>	<b>W9567*</b>		
 <b>1/4 CIRCLE SABRELOC CENTRE POINT SPATULA</b>				9/0 0.3	8/0 0.4	7/0 0.5	6/0 0.7	5/0 1	4/0 1.5	3/0 2
RD-1 	8mm 8mm	30cm	violet				<b>J446G^</b>			
 <b>1/4 CIRCLE ULTIMA SPATULA</b>				10/0 0.2	9/0 0.3	8/0 0.4	7/0 0.5	6/0 0.7	5/0 1	4/0 1.5
CS90-6 	7mm	10cm	violet	<b>V980G^</b>						
 <b>3/8 CIRCLE ULTIMA SPATULA</b>				10/0 0.2	9/0 0.3	8/0 0.4	7/0 0.5	6/0 0.7	5/0 1	4/0 1.5
CS140-6 	6.5mm	10cm	violet	<b>V960G</b>						
 <b>3/8 CIRCLE ULTIMA SPATULA</b>				10/0 0.2	9/0 0.3	8/0 0.4	7/0 0.5	6/0 0.7	5/0 1	4/0 1.5
CS140-6 	6.5mm	10cm	violet		<b>V950G</b>					
CS140-8 	6.75mm 6.75mm	20cm	violet			<b>J547G</b>				
 <b>1/2 CIRCLE ULTIMA SPATULA</b>				10/0 0.2	9/0 0.3	8/0 0.4	7/0 0.5	6/0 0.7	5/0 1	4/0 1.5
CS160-6 	5.5mm 5.5mm	30cm	violet	<b>V448G^</b>						

\* New Zealand Indent

^ Australia & New Zealand Indent

## Coated VICRYL® (Polyglactin 910) Braided Suture

Needle length and description	Length	Colour	Gauge							
			5/0 1	4/0 1.5	3/0 2	2/0 3	0 3.5	1 4	2 5	
<b>NON-NEELED</b>										
LIGAPAK*	135cm 135cm 360cm	violet undyed undyed	<b>J283G*</b>		<b>J205G J285G^</b>	<b>J206G J286G</b>	<b>J207G</b>	<b>J208G</b>		
NON-NEELED	6x45cm	violet					<b>J106T</b>			
NON-NEELED	12x45cm	violet		<b>J103T</b>	<b>J104T</b>	<b>J105T</b>				
NON-NEELED	135cm 135cm 150cm 150cm	violet undyed violet undyed		<b>W9023 J643H*</b>	<b>J614H J606H^</b>	<b>J615H J607H* W9025^</b>	<b>J616H</b>	<b>J617H W9027*</b>	<b>J618H</b>	

\* New Zealand Indent

^ Australia & New Zealand Indent

