


Do your share for the environment

Replace Plastics

The Vernacare Single-Use products are produced using water and recycled over-issued newsprint and the wide range of disposable single use containers are designed to replace re-usable plastic or stainless steel alternatives. Ensuring a clean bedpan, bowl or container every time, proven to help reduce the risk of cross infection and more importantly environmentally friendly.



“It is time for people working in any healthcare situation to think about the validity of single use product.”¹

Vernacare manufacture over **170m** Single-use containers each year

The single-use system is used in over **50** countries worldwide



6 weeks: Vernacare's pulp containers have a similar composition to toilet paper and can fully decompose in just over a month



450+ years: Plastic toileting products have a similar decomposition timeline to plastic bottles making them a significant environmental threat



1. Hospitalhealth.com.au. (2019). Clinical Waste - A Major Environmental Burden.

Vernacare

FOR MORE INFORMATION CONTACT YOUR ACCOUNT MANAGER
OR **FREEPHONE 1800 269 534** OR VISIT www.eboshealthcare.com.au

EBOS
HEALTHCARE

ESSENTIAL TO HEALTH

The Vernacare Single-Use Products



WASHBOWLS

Detergent-proof washbowl can hold warm water and soap for over 2 hours (up to 45 degrees). Unique patented design and formulation.

Washbowl EBOS Code: 118AA100
Case Quantity: 100



MULTI CUP

Ideal for taking urine/stool/sputum samples, dispensing tablets or as a measuring scoop.

EBOS Code: 116AA200
Case Quantity: 200



SMALL BOWL

General purpose bowl, ergonomically shaped with integrated graduations. Ideal to use as vomit bowl.

EBOS Code: 104AA200
Case Quantity: 200



CS RECEIVER

General purpose receptacle, ideal for holding a variety of medical items including dressings, wipes or surgical instruments.

EBOS Code: 106AA300
Case Quantity: 300



KIDNEY DISH

General purpose receptacle, ideal for holding a variety of medical items including dressings, wipes or surgical instruments.

EBOS Code: 105AA260
Case Quantity: 260



MID-STREAM SAMPLER

Ideal for collecting and measuring urine samples.

EBOS Code: 113BA100
Case Quantity: 100



JUG

General purpose jug with ergonomic handle for improved grip and integrated graduations.

EBOS Code: 110AA100
Case Quantity: 100



VERNATRAY

The single use tray makes it easier for clinicians to locate correct syringe and medication.

EBOS Code: 119AA060
Case Quantity: 60

PULP DISPOSAL UNIT: COMPACT

The Vernacare Compact is built for smarter infection prevention. Small enough to be located at the point of care, it provides a quick and hygienic solution for the disposal of Vernacare Single-Use receptacles which have been used to deliver patient care. Inside, the unique SmartFlow™ technology keeps drains flowing.

EBOS Code: 2050ZA001



FEATURING
SMART FLOW™
TECHNOLOGY

Vernacare

FOR MORE INFORMATION CONTACT YOUR ACCOUNT MANAGER
OR **FREEPHONE 1800 269 534** OR VISIT www.eboshealthcare.com.au

EBOS
HEALTHCARE

ESSENTIAL TO HEALTH