



AQ Medical 1000 Vaccine Refrigerator

The AQ Medical 1000 vaccine refrigerator is purpose-built for vaccine storage and follows the Australian Government Department of Health's 'Strive for 5' vaccine storage guidelines.

It was designed with low noise output and consistent temperature readings in mind, and maintains a stable temperature between the recommended +2°C to +8°C. The unit holds safety features such as visual and audible alarms in the event that temperatures fluctuate outside the set range.

All AQ Medical refrigerators are accredited by the Quality Care Pharmacy Program (QCPP), a quality assurance program that was developed by the Pharmacy Guild of Australia in consultation with the Pharmaceutical Society of Australia.

AQ Medical units are also backed by our extensive Avem Quirks service network. With qualified technicians and warehousing facilities across Australia, we can deliver responsive service and exceptional support on a national level.





AQ Medical 1000 Vaccine Refrigerator

Key Features

- QCPP cold chain accredited
- Digital temperature control system preconfigured to hold temperature from +2°C to +8° C in line with national 'Strive for 5' storage guidelines
- Audible and visual auto resetting alarms
- Adjustable shelves
- Lockable door
- Low energy consumption
- Low noise output
- Dedicated technical service line

Product Specifications

Dimensions	Height 1998 mm
	Width 1200 mm
	Depth 697 mm
Gross Weight	kg
Gross Capacity	960 L
Temperature Range	+2°C to +8°C
Number of Shelves	8
Power Consumption/24 hrs.	3.3 kWh
Refrigerant	R600a
Warranty (parts and labour)	24 months

Optional Extras

- Uninterruptible Power Supply unit with back up battery
- AQ Box cloud temperature monitoring system
- Comprehensive programmed maintenance