

# My Practice Essentials

MAY / JUN / JUL 2022

Due to the COVID-19 pandemic please visit [www.eboshealthcare.com.au](http://www.eboshealthcare.com.au) for the fastest way to place your orders.

An EBOS Healthcare Publication  
Primary Care Division Phone: 1800 269 534  
E&OE. Offers not available to medical wholesalers and brokers. All offers valid until: 31st July 2022 unless otherwise specified or while stocks last.  
ALL ENCLOSED PRICES ARE GST EXCLUSIVE

 **EBOS**  
HEALTHCARE

ESSENTIAL TO HEALTH

- 10 % GST will apply to all items unless they are marked with GST FREE
- Freight charges apply for orders under \$250 (GST excluded)
- Prices shown are valid from 1st May 2022 until 31st July 2022

## Your Primary Care Team

- **Business Manager Primary Care** - Moise Wakim  
mwakim@ebosgroup.com.au
- **National Sales Manager** - Matthew Cowlshaw  
mcowlshaw@ebosgroup.com.au

### New South Wales

- **New South Wales Sales Manager** - Carol Wellard  
cwellard@ebosgroup.com.au
- **Key Account Manager**- Ahillen Beadle  
abeadle@ebosgroup.com.au
- **Key Account Manager** - Allison Graham  
agraham@ebosgroup.com.au
- **Key Account Manager** - Bala Nallathambi  
bnallathambi@ebosgroup.com.au
- **Key Account Manager**- Sarah Tiriaki  
stiriaki@ebosgroup.com.au

### Victoria

- **National Corporate Accounts Manager** - Michael Mazengarb  
mmazengarb@ebosgroup.com.au
- **Key Account Manager** - Andrew Brennen  
abrennen@ebosgroup.com.au
- **Key Account Manger** - Nancy Faralla  
nfaralla@ebosgroup.com.au

### Queensland

- **Queensland Sales Manager** - Nick Gower  
ngower@ebosgroup.com.au
- **Key Account Manager** - Ariana Molloy  
amolloy@ebosgroup.com.au
- **Key Account Manager** - Patricio Hernandez  
phernandez@ebosgroup.com.au
- **Key Account Manager** - Shaun Callander  
scallander@ebosgroup.com.au

### Western Australia

- **Western Australia Sales Manager** - Glen Charles  
gcharles@ebosgroup.com.au
- **Key Account Manager** - Deena Doyle  
ddoyle@ebosgroup.com.au
- **Key Account Manager** - Joanna Humphries  
jhumphries@ebosgroup.com.au

### South Australia

- **Key Account Manager** - Janet Pitts  
jpitts@ebosgroup.com.au

## This Issue's Top Picks: Place an order by:

### Examination Couch - 3 Section



### MESI mTABLET System



### BETA 400® F.O. LED Otoscope and Ophthalmoscope Set



Calling us on 1800 269 534



Sending us a Fax on 1800 810 257




Sending us an Email  
customerservice@ebosgroup.com.au




Visiting us online  
www.eboshealthcare.com.au

EBOS Group Australia Pty Ltd  
ABN 38 125 401 247

Availability of products and pricing subject to change without notice.

Free freight on all orders over \$250 (GST Excluded) – excluding equipment indicated with this icon , which will incur additional freight charges. Please call customer service for a quotation.

Prices in catalogue exclude 10% GST. 10% will apply to all items listed unless marked GST FREE. Valid until stocks last and finishes on 31/07/2022.



# Travel is opening up!

See page 3 for **great travel vaccine offers!**

See Pages 3 for details

# Vaccine Essentials:

## FLUCELVAX<sup>®</sup> QUAD

Quadrivalent influenza vaccine (surface antigen, inactivated, prepared in cell cultures)

**Not all influenza vaccines are the same<sup>1</sup>**

Help her make the choice.

### FLUCELVAX<sup>®</sup> QUAD:

- Designed to be an exact match to WHO-selected strains<sup>2-5†</sup>
- †May not predict clinical effects



Recommend FLUCELVAX<sup>®</sup> QUAD for patients 2 years and older<sup>6</sup>

WHO: World Health Organisation

**PBS Information:** This product is not listed on the National Immunisation Program (NIP) or the PBS.

Before prescribing, please review Product Information available at [www.seqirus.com.au/products](http://www.seqirus.com.au/products)

**MINIMUM PRODUCT INFORMATION Fluclavax<sup>®</sup> Quad** (Inactivated Quadrivalent Influenza Vaccine (prepared in cell cultures)) - Suspension for injection. **Indication:** For the prevention of influenza caused by Influenza Virus, Types A and B contained in the vaccine, for use in adults and children 2 years of age and older. **Contraindications:** Individuals with known severe allergic reactions (e.g. anaphylaxis) to any component of the vaccine or to a previous dose of any influenza vaccine. **Precautions:** Treatment and supervision for anaphylactic reactions should be available. Antibody response in immunocompromised patients may be lower. Careful consideration for use in history of Guillain-Barré syndrome within 6 weeks of previous influenza vaccination. Postpone immunisation in patients with febrile illness or acute infection. A protective immune response may not be elicited in all vaccine recipients. **Use in Pregnancy (Category B1):** The safety of Fluclavax<sup>®</sup> Quad in pregnancy has not been assessed in clinical trials. There are no reproductive and developmental toxicity studies with Fluclavax<sup>®</sup> Quad. Reproductive and developmental toxicity data from cell-based trivalent influenza vaccine (TIVc) do not predict an increased risk of developmental abnormalities. **Use in Lactation:** Not been evaluated in nursing mothers. **Interactions:** No data available on co-administration with other vaccines. Based on clinical experience with TIVc, Fluclavax<sup>®</sup> Quad can be given at the same time as other vaccines (in separate limbs) (adverse reactions may be intensified). **Adverse Effects:** Local - pain, erythema. Systemic - headache, fatigue, irritability (children 2-6 years), upper respiratory tract infection, nasopharyngitis (See full PI). **Dosage and Administration:** Intramuscular injection only: single 0.5mL dose (children <9 years may require 2 doses - see full PI). **Presentation:** 0.5mL single-use syringe. Storage: Store at 2°C to 8°C; must not be frozen; protect from light. Based on Approved Product Information. **References:** 1. ATAGI. Statement on the administration of seasonal influenza vaccines in 2022. [www.health.gov.au/resources/publications/atagi-advice-on-seasonal-influenza-vaccines-in-2022](http://www.health.gov.au/resources/publications/atagi-advice-on-seasonal-influenza-vaccines-in-2022) Accessed March 2022. 2. Rajaram S *et al.* *Ther Adv Vaccines Immunother.* 2020;8:2515135520908121. 3. CDC. Cell-based flu vaccines. Available from: [www.cdc.gov/flu/prevent/cell-based.htm](http://www.cdc.gov/flu/prevent/cell-based.htm) Accessed March 2022. 4. Mabrouk T *et al.* *Dev Biol (Basel).* 2002;110:125-134. 5. WHO. Influenza A(H3N2) lineage cell culture-derived CVVs for development and production of vaccines for use in the 2022 SH influenza season; Available from: [https://cdn.who.int/media/docs/default-source/influenza/cvvs/cvv-southern-hemisphere-2022/a\(h3n2\)--cell-culture-derived.pdf](https://cdn.who.int/media/docs/default-source/influenza/cvvs/cvv-southern-hemisphere-2022/a(h3n2)--cell-culture-derived.pdf) Accessed March 2022. 6. FLUCELVAX<sup>®</sup> QUAD Approved Product Information. Seqirus (Australia) Pty Ltd Melbourne, Victoria. ABN: 66 120 398 067. [www.seqirus.com.au](http://www.seqirus.com.au). Seqirus Medical Information: 1800 642 865. FLUCELVAX<sup>®</sup> is a registered trademark of Seqirus UK or its affiliates. Seqirus<sup>™</sup> is a trademark of Seqirus UK Limited or its affiliates. Date of preparation: March 2022. ANZ-QIVc-22-0027. 001351-A.

**Seqirus**<sup>™</sup>  
A CSL COMPANY

<b>Trumenba<sup>®</sup></b> Meningococcal Group B Vaccine PF33272	\$147 <b>5+</b> \$140	<b>Priorix<sup>®</sup></b> Measles, mumps, rubella vaccine, live SKLN1359	\$27.90 <b>10+</b> \$22
<b>BEXSERO<sup>®</sup></b> Meningococcal Group B Vaccine SKN1010	\$115 <b>10+</b> \$97	<b>ViVAXIM<sup>®</sup></b> Combined Purified V Polysaccharide Typhoid and Inactivated Hepatitis A Vaccine AV8090	\$143 <b>10+</b> \$133
<b>STAMARIL<sup>®</sup></b> Yellow Fever Vaccine (Live, Stabilised) AV8040	\$103 <b>10+</b> \$88	<b>VAQTA<sup>®</sup></b> Hepatitis A vaccine, inactivated CS80600001	\$49 <b>10+</b> \$42
<b>Rabipur<sup>®</sup></b> (Inactivated rabies virus vaccine) CS08510101	\$119 <b>5+</b> \$102	<b>Adacel<sup>®</sup></b> Pentavalent Acellular Diphtheria, Tetanus, Pertussis, Polysaccharide Acellular Pertussis and Polysaccharide Hib Vaccine AV350693	\$35 <b>5+</b> \$28
<b>TYPHIM Vi<sup>™</sup></b> Salmonella Typhi Vi Polysaccharide Vaccine AV8020	\$58 <b>10+</b> \$52	<b>DUKORAL<sup>®</sup></b> Cholera Heat Stable Bacterin-vaccine AV8170	Single \$72 <b>10+</b> \$57
<b>I POL<sup>®</sup></b> Inactivated polio vaccine (IIP) CS08160101	\$44 <b>10+</b> \$37	<b>DUKORAL<sup>®</sup></b> Cholera Heat Stable Bacterin-vaccine AV8172	Twin Pack \$145 <b>10+</b> \$114
<b>M-M-R II<sup>®</sup></b> Measles, Mumps, & Rubella Virus Vaccine (Live) CS07790101	\$39 <b>10+</b> \$31	<b>havrix<sup>®</sup></b> Inactivated hepatitis A virus vaccine SKN3708	Adult \$53 <b>10+</b> \$39
<b>Prevenar 13<sup>®</sup></b> Pneumococcal polysaccharide vaccine AY01461	\$149 <b>10+</b> \$119	<b>havrix<sup>®</sup></b> Inactivated hepatitis A virus vaccine SKN3718	Junior \$39 <b>10+</b> \$31
<b>AVAXIM<sup>™</sup></b> Inactivated Hepatitis A Vaccine AV8050	\$55 <b>10+</b> \$48	<b>ADACEL<sup>®</sup> polio</b> AV8190	\$62 <b>10+</b> \$56
<b>boostrix<sup>®</sup></b> Combined diphtheria-tetanus-acellular pertussis (dTpa) vaccine SKN1014	\$35 <b>10+</b> \$29	<b>MENVEO<sup>®</sup></b> Meningococcal Groups A, C, Y and W-135 Polysaccharide CRM PP Conjugate Vaccine SKN4300	\$79 <b>10+</b> \$59
<b>Menactra<sup>®</sup></b> Meningococcal (Groups A, C, Y and W-135) Polysaccharide Diphtheria Toxoid Conjugate Vaccine FI324333	\$80 <b>10+</b> \$59	<b>twinrix<sup>®</sup></b> combined hepatitis A and hepatitis B vaccine SKN3186	Adult \$68 <b>10+</b> \$58
<b>Nimenrix<sup>®</sup></b> Meningococcal polysaccharide groups A, C, Y and W-135 conjugate vaccine SKN3151	\$62 <b>10+</b> \$54	<b>twinrix<sup>®</sup></b> combined hepatitis A and hepatitis B vaccine SKN3175	Junior \$48 <b>10+</b> \$39



- 10 % GST will apply to all items unless they are marked with GST FREE
- Freight charges apply for orders under \$250 (GST excluded)
- Prices shown are valid from 1st May 2022 until 31st July 2022

## Refrigeration Essentials:

EXCEPTIONAL QUALITY AND TRULY OUTSTANDING FEATURES, AT AN AFFORDABLE PRICE.

### Key Features:

- QCPP Compliant
- Independently tested by The Pharmacy Guild of NSW
- Touch Screen Digital Controller
- Audible & Visual Alarm
- Built In Lock
- High Ambient Rating 38°C
- 24 months, 60 months full warranty if registered online and serviced and calibrated annually



**\$3,595**

EBOS Code: VACCSAFE4



**\$5,799**

EBOS Code: VACCSAFE6



**\$7,395**

EBOS Code: VACCSAFE10

Premium  
**Vacc-Safe®**

### Specifications:

Model	VACCSAFE4	VACCSAFE6	VACCSAFE10
<b>Dimensions External (W x D x H) mm</b>	606 x 620 x 1940	700 x 790 x 2058	1130 x 720 x 2023
<b>Gross Litres</b>	400L	600L	1000L
<b>Compliance</b>	QCPP and Cold Chain compliant + Pharmacy Guild Tested	QCPP and Cold Chain compliant + Pharmacy Guild Tested	QCPP and Cold Chain compliant + Pharmacy Guild Tested
<b>Warranty</b>	24 months, 60 months full warranty if registered online & serviced & calibrated annually	24 months, 60 months full warranty if registered online & serviced & calibrated annually	24 months, 60 months full warranty if registered online & serviced & calibrated annually
<b>Weight</b>	75Kg	110Kg	130Kg
<b>Operating range</b>	2.C to 8.C	2.C to 8.C	2.C to 8.C
<b>Shelves</b>	7	6	12 (6 per door)
<b>Lock</b>	Self closing lockable doors	Self closing lockable doors	Self closing lockable door. Door lock in top panel
<b>Alarm</b>	NO/NC External Alarm Contacts Fitted	NO/NC External Alarm Contacts Fitted	NO/NC External Alarm Contacts Fitted
<b>Power</b>	10 AMP	10 AMP	10 AMP
<b>Monitoring</b>	Medi Log 2 Data Logger Audible and Visible Alerts	Medi Log 2 Data Logger Audible and Visible Alerts	Medi Log 2 Data Logger Audible and Visible Alerts
<b>Cooling</b>	Fan Assisted	Fan Assisted	Fan Assisted
<b>Body Type</b>	White Powder Coated Steel Internal and External	White Powder Coated Steel Internal and External	White Powder Coated Steel Internal and External
<b>Power Consumption</b>	4.935KWH (Over 24hrs)	6.209KWH (Over 24hrs)	9.903KWH (Over 24hrs)

# Refrigeration Essentials:

COMBINING SUPERB PERFORMANCE AND QUALITY WITH GREAT VALUE.

## Key Features:

- QCPP Compliant
- Independently tested by the Pharmacy Guild Of NSW
- Self-closing glass door with inbuilt lock
- 24-month comprehensive warranty
- Rated to 38 Degrees



**\$2,850**

EBOS Code: VACCMED4



**\$4,495**

EBOS Code: VACCMED6



**\$6,895**

EBOS Code: VACCMED10

## Specifications:

Model	VACCMED4	VACCMED6	VACCMED10
<b>Dimensions External (W x D x H) mm</b>	606 x 620 x 1940	700 x 790 x 2058	11130 x 720 x 2023
<b>Gross Litres</b>	400L	600L	1000L
<b>Compliance</b>	QCPP and Cold Chain compliant + Pharmacy Guild Tested	QCPP and Cold Chain compliant + Pharmacy Guild Tested	QCPP and Cold Chain compliant + Pharmacy Guild Tested
<b>Warranty</b>	24 months parts and labour	24 months parts and labour	24 months parts and labour
<b>Weight</b>	75Kg	110Kg	142Kg
<b>Operating range</b>	2.C to 8.C	2.C to 8.C	2.C to 8.C
<b>Shelves</b>	7	6	12
<b>Lock</b>	Self closing lockable doors	Self closing lockable doors	Self closing lockable doors
<b>Alarm</b>	High / Low and Door Open Audible and Visual Alarm	High / Low and Door Open Audible and Visual Alarm	High / Low and Door Open Audible and Visual Alarm
<b>Power</b>	10 AMP	10 AMP	10 AMP
<b>Monitoring</b>	Medi Log 2 Data Logger Audible and Visible Alerts	Medi Log 2 Data Logger Audible and Visible Alerts	Medi Log 2 Data Logger Audible and Visible Alerts
<b>Cooling</b>	Fan Assisted	Fan Assisted	Fan Assisted
<b>Body Type</b>	White Powder Coated Steel Internal and External	White Powder Coated Steel Internal and External	White Powder Coated Steel Internal and External
<b>Power Consumption</b>	4.935KWH (Over 24hrs)	6.209KWH (Over 24hrs)	9.903KWH (Over 24hrs)

- 10 % GST will apply to all items unless they are marked with GST FREE
- Freight charges apply for orders under \$250 (GST excluded)
- Prices shown are valid from 1st May 2022 until 31st July 2022

## Cryosurgery Essentials:

### CryoPro Maxi Flask 500mL

The CryoPro Cryosurgical Flask is robust in design and construction, providing a reliable tool for your surgery. It builds in safety features like a relief valve and a precision trigger for ease of use.

Flask comes with a set of 6-tips and probes



**\$875**

EBOS Code: ELCPRMAX

### CryoPro Mini Flask 350mL

- Suitable for use by both GP and Specialists
- 350ml Capacity
- Comes with applicator set (5 tips and 1 cleaning adaptor)
- 3 year guarantee
- Parts can be universal and interchangeable with Brymill



**\$830**

EBOS Code: ELCPRO

### Dewar Liquid Nitrogen Container XT 20 Litre

- Superior vacuum performance with super insulation
- Evaporation rate and static holding time are nominal.
- Work time is an arbitrary, reference-only value to estimate container performance under the actual operating conditions.
- Comes with 5 year vacuum warranty, 2 year container warranty
- Each unit comes complete with safety lock lid, neck, plug, leather (PU) bag and ladle



RRP: \$ 1,295

**\$959**

SAVE  
**\$336**

EBOS Code: TW20XT0

### Ladle for Dewar XT with handle



**\$99**

EBOS Code: TWDIP3

### Polar Grip™ Size 9 Pair



**\$29.95**

EBOS Code: AN23700

## Blood Pressure Essentials:

### Ambulatory Blood Pressure Monitor TM-2441

- Light weight – 122 grams (excluding batteries)
- Oscillometric measuring method
- AFib+ Irregular Heartbeat detection (IHB) and AFib screening
- Multiple measurement intervals – adapts to patient's and healthcare professional's exact requirements
- 3 years warranty



**Item number 11607**  
**Fee: \$107.20**  
**Benefit:**  
**75% = \$80.40**  
**85% = \$91.15**

EBOS Code	Description	UOM	Box/50
200693	SMALL CUFF (LEFT) TM-2440-1	Each	<b>\$155</b>
200690	ADULT CUFF (LEFT) TM-2440-1	Each	<b>\$165</b>
200692	LARGE CUFF (LEFT) TM-2440-1	Each	<b>\$210</b>
200694	XLARGE CUFF (LEFT) TM-2440-1	Each	<b>\$230</b>

**\$1,895**

EBOS Code: 200695

### Ambulatory Blood Pressure Monitor TM-2440

- Automatic or manual measurement method
- Oscillometric measuring
- NEW pulse correction algorithm for faster measurement & noise rejection
- Multiple measurement intervals – adapts to patient requirements
- A/A accuracy, BHS protocol validated
- Includes 20-31cm cuff



**Item number 11607**  
**Fee: \$107.20**  
**Benefit:**  
**75% = \$80.40**  
**85% = \$91.15**

EBOS Code	Description	UOM	Box/50
200693	SMALL CUFF (LEFT) TM-2440-1	Each	<b>\$155</b>
200690	ADULT CUFF (LEFT) TM-2440-1	Each	<b>\$165</b>
200692	LARGE CUFF (LEFT) TM-2440-1	Each	<b>\$210</b>
200694	XLARGE CUFF (LEFT) TM-2440-1	Each	<b>\$230</b>

**\$1,495**

EBOS Code: ANDTM2440

- 10 % GST will apply to all items unless they are marked with GST FREE
- Freight charges apply for orders under \$250 (GST excluded)
- Prices shown are valid from 1st May 2022 until 31st July 2022

## HEINE Essentials:

### HEINE DELTAone

Mobility and quality have now become one.  
The HEINE DELTAone Dermatoscope.



This mobile masterpiece of dermatoscopy is taking medical professionals by storm. The crystal-clear image was accomplished using the new achromatic HEINE optical system. What's more, our LEDHQ produces extraordinary colour rendering. Plus the polarisation provides glare-free and non-reflective working conditions. The new DELTAone is a tool you'll gladly carry in your pocket.



**21mm  
Field of  
View**

EBOS Code	Description	UOM	RRP	PRICE	SAVE
HEK2102830	DELTAone Standard Edition	Each	\$1,969	<b>\$1,799</b>	<b>\$170</b>
HEK2112830	DELTAone White Edition	Each	\$1,969	<b>\$1,799</b>	<b>\$170</b>

### HEINE DELTA 30

If there's something to discover, you'll be sure to see it. The HEINE DELTA 30 Dermatoscope.

You'll see things you've never seen before. Our new high-end dermatoscope provides the best imaging we've ever developed. Combined with a field of view measuring a veritable 30 millimetres and unprecedented colour rendering, in LEDHQ of course. Polarisation provides glare-free and non-reflective working conditions. Head off on a journey of discovery.

**30mm  
Field of  
View**



EBOS Code	Description	UOM	RRP	PRICE	SAVE
HEK2302830	DELTA 30 Dermoscope	Each	\$2,699	<b>\$2,499</b>	<b>\$200</b>
HEX0029921	DELTA 30 Table Charging Station	Each	\$429	<b>\$399</b>	<b>\$30</b>

### HEINE GAMMA G5®

The HEINE standard sphygmomanometer. Virtually the same as the G7 Model except the housing is moulded of 2 component thermo plastic. 5 year manufacturers warranty



**5  
YEAR  
GUARANTEE**

RRP: \$186

**\$135** **SAVE \$51**

EBOS Code: HEM0009235





## HEINE Essentials:

### BETA 400® F.O. LED Otoscope & Ophthalmoscope Set

- LED in HEINE Quality in fibre optic illumination for the highest colour fidelity and edge to edge illumination
- Exclusive 4.2x magnification
- Multi coated precision optics
- AOS technology
- Patented, unique, stepless dimming from 3 % to 100 %
- Range of lenses: 27 lenses from 35D to +40D
- No stray light due to the recessed, multi coated viewing window
- Ergonomic design fits the orbita comfortably, and blocks out ambient light
- Charging options: Desktop, USB and Battery

**Set includes:**

- 1 set of reusable tips
- 10 Disposable tips
- Hard case
- BETA 4 NT rechargeable handle
- NT 4 table charger
- 2nd handle included FREE



RRP: \$3,156

**\$2,299**

SAVE  
**\$857**



EBOS Code: HEA1532442

### K180® Diagnostic Set

- Large field of view
- K180 Ophthalmoscope
- K180 F.O. Otoscope
- 1 set of reusable tips
- Auto on-off recharge feature



RRP: \$1,113

**\$825**

SAVE  
**\$288**



EBOS Code: HEA2791011

### HEINE mini 3000® Diagnostic Set

- No stray light due to the recessed, multi-coated viewing window
- mini 3000 Ophthalmoscope
- mini 3000 Otoscope
- 2 mini 3000 battery handles
- 1 set of reusable tips
- 10 disposable ear tips
- Hard case



RRP: \$729

**\$549**

SAVE  
**\$180**



EBOS Code: HED8732102

- 10 % GST will apply to all items unless they are marked with GST FREE
- Freight charges apply for orders under \$250 (GST excluded)
- Prices shown are valid from 1st May 2022 until 31st July 2022

## Furniture Essentials: Emergency Trolley

- Defibrillator shelf
- Extension table
- 304 Stainless steel IV pole
- Cardiac CPR board
- Waste bin
- External power holder and plug board
- Oxygen tank holder
- 304 Stainless steel wire basket
- Central or seals lock
- Castors with brakes
- Non slip mat
- Luxurious push handle
- Soft closing drawers
- Dimensions: 650 x 480 x 980mm
- Drawer Sizes: 3 x 110mm, 1 x 135mm, 1 x 195mm
- G.W 55kg
- Central locking system through side key



RRP: \$1,795

**\$1,384** **SAVE \$411**

EBOS Code: PMEMTR



**All  
Accessories  
included!**

## Free Standing Couch

- Fixed height 82cm
- New easy lift and lock mechanism for easy assembly
- Lift up back section for assessment
- Comes flat packed
- Total Length: 195cm
- Width: 71cm
- Height: 82cm
- Loading weight: 280Kgs
- Black and Navy Blue also available



RRP: \$487

**\$345** **SAVE \$142**

EBOS Code: PMFSGR



## Powder Coated Trolley's

- Bottom shelf
- Soft closing magnetic drawer
- Smooth running lockable castors
- Top and bottom side rails on three sides
- 125mm deep drawer
- 95mm front castors
- Dimensions: 50 x 50 x 96.7cm (LxWxH)
- Powder coated finish

### 1 Drawer

- 500 x 500 x 900mm



RRP: \$365

**\$309** **SAVE \$56**

EBOS Code: PMPCT01



### 2 Drawer

- 500 x 500 x 900mm



RRP: \$375

**\$359** **SAVE \$16**

EBOS Code: PMPCT02



### 3 Drawer

- 500 x 500 x 900mm



RRP: \$445

**\$409** **SAVE \$36**

EBOS Code: PMPCT03



## Furniture Essentials: Examination Couch - 3 Section

- Electric Hi Lo operation by foot switch
- Head section face hole and plug
- Dual gas strut supported head and back sections -45° to +70°
- Retractable castors for safety and mobility
- 5 Year frame and 3 year electrics warranty

### Specifications:

Length	
Total	195cm
Head Section	63cm
Middle Section	63cm
Leg Section	63cm
Width	
	71cm
Height Adjustable	46.5cm - 92.5cm
Loading Weight	225kg
Padding thickness	5cm
Foam Density	36kg/m3

**BONUS  
Round  
Stool!**



Scan to see the  
examination couch  
in action.

**\$1,259**

EBOS Code: ET33NB



## Stainless Steel Trolley's

- High grade 304 stainless steel throughout
- Soft closing magnetic drawers
- German made ball bearing drawer runners
- 125mm deep drawers
- 95 mm smooth rolling castors for easy movability
- Locks on front castors
- 125mm deep drawers
- Top and bottom side rails on all three sides
- Double safety folds
- Smooth welding

### 1 Drawer

- 500 x 500 x 900mm



RRP: \$ 495

**\$409** SAVE \$86

EBOS Code: PMSST01



### 2 Drawer

- 500 x 500 x 900mm



RRP: \$695

**\$515** SAVE \$180

EBOS Code: PMSST02



### 3 Drawer

- 500 x 500 x 900mm



RRP: \$775

**\$669** SAVE \$106

EBOS Code: PMSST03



- 10 % GST will apply to all items unless they are marked with GST FREE
- Freight charges apply for orders under \$250 (GST excluded)
- Prices shown are valid from 1st May 2022 until 31st July 2022

## MESI Essentials:

### MESI mTABLET Spirometry Module

- Multiple measurement modes and parameter calculation options
- Automatic selection of the best of all repeated maneuvers with the BestBreath™ detection
- Detailed report with the ability to switch between charts and values for a clear interpretation
- Option to compare different measurements on the same screen through MESI mRECORDS
- Pneumatograph technology with integrated self-calibration for accurate measurements at any time

**\$2,195**

EBOS Code: SPIROSYSCS

#### Accessories

EBOS Code	Description	UOM	PRICE
SPIROMDMP	MESI DISPOSABLE FLOW TRANSDUCER	Box/100	<b>\$285</b>
SPIROMDFIL	MESI MTABLET SPIRO FILTER	Box/50	<b>\$125</b>
MTABTRLLY	MESI MTABLET TROLLEY	Each	<b>\$1,700</b>



\*MESI mTABLET not included

**MESI**  
Simplifying Diagnostics



## MESI Essentials:



Simplifying Diagnostics

Provide quality healthcare from a safe distance  
**with the MESI mTablet System**



Product Code	Description	Unit	Quantity	Unit Price (excl GST)	Package Deal Price
MTABSYSABI	MESI MTABLET ABI DIAG SYSTEM	each	1	\$4,800.00	\$4,800.00
ECGSYS	MESI MTABLET 12-LEAD ECG DIAG MODULE	each	1	\$2,900.00	\$2,900.00
MTABABICL	MESI MTABLET ABI CUFF SET (4) LARGE	each	1	\$450.00	FREE
AVBP	AVERAGE BLOOD PRESSURE APP	each	1	\$400.00	FREE
PHOTO	PHOTO APP	each	1	\$350.00	FREE
BPAPP	BLOOD PRESSURE APP	each	1	\$400.00	FREE
ININSTALL	MESI FULL INTEGRATION INSTALL	each	1	\$400.00	FREE
<b>TOTAL</b>				<b>\$9,700.00</b>	<b>\$7,700.00</b>

- 10% GST will apply to all items unless they are marked with GST FREE
- Freight charges apply for orders under \$250 (GST excluded)
- Prices shown are valid from 1st May 2022 until 31st July 2022

## Practice Essentials:

### Ordisi Mobile Examination Light

Cool light source

- Expected LED lifetime 50,000 hours
- Very low power consumption (10w)
- Strong, white light
- 42,000 lux/0.5m
- 12,000 lux/1m
- Lightweight mobile base
- 12 month warranty



**EOFY Special!**

RRP: \$ 758

**\$499**

SAVE \$259



EBOS Code: ODLEDMEXAM

### Ordisi Wall Examination Light

Includes Swivel Bracket

- Cool light source
- Expected LED lifetime 50,000 hours
- Very low power consumption (10w)
- Strong, white light
- 42,000 lux/0.5m
- 12,000 lux/1m
- Goose style neck for easy maneuverability
- 12 month warranty



**EOFY Special!**

RRP: \$ 763

**\$479**

SAVE \$284



EBOS Code: ODFLH2731W

### Projet 101 Pulsatile Ear Irrigator

Designed by Doctors with innovative pulsatile irrigation



**\$419**

SAVE \$26

2 year warranty

EBOS Code: PROJET101

EBOS Code	Description	UOM	PRICE
PROJET-TIP	PROJET TIPS	Box/100	\$142

### Seca 354 Electronic Baby Scales

- Scale Capacity: 20kg
- Scale Graduation: 10g < 10kg > 20g
- Functions: TARE, auto-HOLD, BMIF, kg/lbs
- switch-over



**\$185**

EBOS Code: PN354



### seca 813

- Capacity: 200kg
- Graduation: 100g
- Functions: auto-HOLD, kg/lb switch over, tip-on automatic switch-on, automatic switch-off.



**\$175**

EBOS Code: PN813



### Retractable Measuring Tape 1.5m



**\$2.95**


EBOS Code: SAS423

## Wound Care Essentials:


HARTMANN



**simply**  
wound care



Scan the QR code with your phone  
or search for Simply Wound Care



- 10 % GST will apply to all items unless they are marked with GST FREE
- Freight charges apply for orders under \$250 (GST excluded)
- Prices shown are valid from 1st May 2022 until 31st July 2022

## Hand Hygiene Essentials



**Whiteley**   
Protecting lives and livelihoods since 1933

### Preventing infection transmission shouldn't mean dry skin.

Control infectious disease whilst preventing skin moisture loss.

- ✓ Leaves hands hygienic, soft and replenished
- ✓ Contains natural emollients
- ✓ Maintains skin integrity

Bactol® Clear is an antibacterial hand rub for use without water. It is colour free and pH balanced, leaving skin feeling soft and refreshed.

Proudly Australian made, you can rely on Bactol® Clear to protect your patients, your staff and your reputation.

**Always read the label. Follow the directions for use.**

**Bactol® Clear available in:**

- ✓ 500mL
- ✓ 1 Litre pods for use in Whiteley dispensing systems.



AUST R: 335834

For more information call **1800 833 566** or go to [whiteley.com.au](http://whiteley.com.au)

© Registered Trademark Whiteley Corporation Pty. Ltd. © 2022

## Bactol® Clear Antibacterial Hand Rub

Bactol® Clear is colour free antibacterial hand rub that kills bacteria without water. It maintains skin integrity with a pH balanced formula including natural emollients that prevent moisture loss leaving your hands hygienic, soft and replenished. As recommended by WHO, Bactol Clear is formulated with an alcohol concentration of 70% Ethanol v/v, for maximum efficacy.



EBOS Code	Description	UOM	PRICE
WH012132	BACTOL CLEAR 500mL	Each	<b>\$14.10</b>
WH012134	BACTOL CLEAR 1000mL	Box/6	<b>\$157.00</b>
WH302017	BACTOL HAND GEL 1L DISPENSER AUTO	Each	<b>\$92.00</b>
WH302016	BACTOL HAND GEL 1L DISPENSER MANUAL	Each	<b>\$39.00</b>



## Woundcare Essentials:

**3M** Science.  
Applied to Life.™

 Non-adhesive



**3M™ Coban™ NL Non-Latex  
Self-Adherent Wrap**

Other non-adhesive 3M medical  
securement solutions: 3M™ Coban™  
Self-Adherent Wrap (contains latex)



**3M™ Medical Tapes and Wraps  
Transforming Patient Care**

 Multi-purpose



**3M™ Micropore™  
Surgical Tape**     **3M™ Micropore™  
S Surgical Tape**

Other multi-purpose 3M medical  
securement solutions: 3M™ Transpore™  
Surgical Tape

## 3M Micropore Surgical Tape

3M Micropore Surgical Tape is a gentle, general purpose paper tape which is highly breathable to maintain skin integrity.

- Latex-free and hypoallergenic for sensitive patients
- Highly breathable to maintain skin integrity
- Holds well on damp skin for secure placement

EBOS Code	Description	UOM	PRICE
3M1530-0	12mmx9.1m	Box/24	<b>\$12.00</b>
3M1530-1	25mmx9.1m	Box/12	
3M1530-2	50mmx9.1m	Box/6	
3M1530-3	72mmx9.1m	Box/4	



 Flexible



**3M™ Medipore™ H  
Soft Cloth Surgical Tape**

Other flexible 3M medical  
securement solutions:  
3M™ Medipore™ H Surgical Tape

 High-strength



**3M™ Durapore™  
Surgical Tape**

Individually packaged single-patient-  
use rolls: Can help protect patients by  
reducing the risk of cross-contamination



## Medipore H Soft Cloth Surgical Tape

- Gentle, breathable to maintain skin integrity
- Easy-tear perforated rolls
- Excellent cross and diagonal stretch to accommodate swelling and movement
- Hypoallergenic and not made with natural rubber latex
- Water resistant

EBOS Code	Description	UOM	PRICE
3M2961	25mmx9.1m	Box/24	<b>\$129.30</b>
3M2962	50mmx9.1m	Each	<b>\$10.80</b>
3M2963	75mmx9.1m	Box/12	<b>\$193.90</b>



Medical Solutions Division. Available to Healthcare Professionals only.  
3M Australia Pty Limited. ABN 90 00 100 096. Building A, 1 Rivett Road, North Ryde, NSW 2113. Ph: 1300 363 454 3M.com.au/medical. 3M, "3M Science. Applied to Life.", 3M, Coban, Transpore, Medipore, Durapore, Micropore and Microfoam are trademarks of 3M Company. 3M 2022. All rights reserved.

- 10 % GST will apply to all items unless they are marked with GST FREE
- Freight charges apply for orders under \$250 (GST excluded)
- Prices shown are valid from 1st May 2022 until 31st July 2022

## inhealth™ Essentials:

### inhealth™ Examination Gloves Nitrile

- Ambidextrous
- Textured fingertips
- Tested for use with chemotherapy drugs



EBOS Code	Description	Box/100
EBNGS150	SMALL	<b>\$9.90</b>
EBNGM150	MEDIUM	
EBNGL150	LARGE	
EBNGXL150	XLARGE	



### inhealth™ 4ply Underpads

The underpads provides soft, non-woven fabric with a highly absorbent fill material to maximise ability to absorb unwanted fluids. The underpads contains an absorbent inner to lock away wetness. It is sealed on four sides, with a waterproof backing to prevent strikethrough.

- Absorbent – Absorbs and holds in liquid.
- Soft top cover fabric to increase comfort while maintaining its strength
- Soft waterproof backing is silent and minimise movement when in use.
- Waterproof backing prevents leaking



### inhealth™ IV Pressure Pad Haemostatic



**\$27.50**

EBOS Code: CCIV36



### inhealth™ Bed Roll

Inhealth bed rolls are 49cm wide x 41.5 meters in length with perforations every 42cm. Within a clinical situation, the Bed Roll is predominantly used in the general practice examination room, however other uses include, surface protector, surface wipe, non-sterile procedural area towel and general cleaning.

UOM: Each



**\$10.50**

EBOS Code: KLCRL049



EBOS Code	UOM	Pkt/300
CCUP4060	Box/200	<b>\$42.50</b>



## inhealth™ Essentials:

### Interleaved Hand Towel

- Fits most dispensers



EBOS Code	Description	UOM	PRICE
KLITD19	90 Shts 29.5 x 19cm	Ctn/24	<b>\$42.50</b>

### All Purpose Towel - Low Lint



EBOS Code	Description	UOM	PRICE
CCAPT3033	All purpose towel	Box/75	<b>\$6.50</b>

### Carbon Steel Scalpel Blades



EBOS Code	Description	UOM	PRICE
PCS10	Size 10	Box/100	<b>\$21.40</b>
PCS11	Size 11	Box/100	
PCS15	Size 15	Box/100	
PCS20	Size 20	Box/100	
PCS22	Size 22	Box/100	
PCS23	Size 23	Box/100	



### Disposable Scalpel Blades and Handles



EBOS Code	Description	UOM	PRICE
PDSS10	Size 10	Box/10	<b>\$11.90</b>
PDSS11	Size 11	Box/10	
PDSS15	Size 15	Box/10	
PDSS23	Size 23	Box/10	



### Carbon Steel Stitch Cutters



EBOS Code	Description	UOM	PRICE
PSSCS	Regular	Box/100	<b>\$23.90</b>



- 10 % GST will apply to all items unless they are marked with GST FREE
- Freight charges apply for orders under \$250 (GST excluded)
- Prices shown are valid from 1st May 2022 until 31st July 2022

## Med Con Essentials:

# Hypafix®



- Hypafix® is a self-adhesive, non-woven fabric sheet for secure and rapid fixation.
- It is made from white polyester material and coated with hypoallergenic adhesive on quick-release backing paper.
- Hypafix also provides complete coverage of the entire dressing, thus reducing the risk of contamination.
- Its excellent conformability simplifies shaping around joints and awkward body contours to allow the patient greater freedom of movement.

EBOS Code	Description	EACH
SN4208	HYPAFIX 2.5cmx10m	<b>\$9.55</b>
SN4209	HYPAFIX 5cmx10m	<b>\$13.80</b>
SN4210	HYPAFIX 10cmx10m	<b>\$23.00</b>
SN4211	HYPAFIX 15cmx10m	<b>\$34.95</b>



## Instant-view® One-Step Pregnancy Test

- 10 mIU/ml sensitivity
- Simple and easy to use one-step testing
- Rapid results may be seen as quickly as 1-3 minutes



**\$24.00**

EBOS Code: NIAPT25

## TUFFIE 5

### Cleaning and Disinfecting Wipes

Aids the removal of organic matter and improves infection prevention and control.



EBOS Code	Description	UOM	PRICE
901SW225	Universal Wipes	Tub/225	<b>\$19.95</b>



## TUFFIE

### Alcohol Free Detergent Wipes

- Improved cleaning performance
- Non-streak formula
- Faster drying time
- Active Odour Control



EBOS Code	Description	UOM	PRICE
901DW225	Detergent Wipes	Tub/225	<b>\$18.95</b>

