

# PR PHARMACY REFRIGERATOR



## FEATURES

- Fan Forced System
- Top Mounted Door Lock
- Interior LED
- Double Glazed Door
- Hi/Lo Readout
- Door Fan Switch
- Digital Display with Built-in Alarms
- Highly Energy Efficient Refrigeration System
- QCPP Compliant
- ARTG No. 478745



# PR PHARMACY REFRIGERATOR



## SPECIFICATIONS

Model	PR150	PR350
Gross Capacity	135 L	350 L
External Dimensions in mm (w/d/h including castors)	595 x 650 x 830	595 x 650 x 1850
Internal Dimensions in mm (w/d/h)	455 x 445 x 665	455 x 445 x 1685
Packed Dimensions in mm (w/d/h including pallet)	640 x 670 x 1005	640 x 670 x 2010
Gross Weight - kg (packed weight)	51	94
Net Weight - kg	45	77
Energy Consumption	0.86 kWh/24	1.1 kWh/24
Voltage	240V	
Temperature Range	2°C to 8°C	
Ambient Temperature Range	10°C to 32°C	
Refrigerant	R600a	
Cooling System / Defrost	Fan Forced / Automatic	
Exterior Cabinet Finish	Enamelled Steel	
Inner Liner Material	ABS	
Door Hinge	RH	
Type of Controller	Digital	
Temperature Display	External	
Door Ajar Alarms	Standard	
Fan Cutout on Opening	Standard	
High and Low Readout	Visual / Audible	
Power Requirement	10amp GPO	
Amps	0.6	0.8
Castors	Levelling Legs	
Adjustable Shelves	4	7
Usable Shelf Area in mm (w/d)	3 x 530 x 652 1 x 438 x 202	6 x 455 x 405 1 x 438 x 202
Light	LED	
Self Closing Door	Standard on RH Only	
Data Logger	Optional	
Uninterruptible Power Supply Unit with Back Up Battery	Optional	
Preventative Maintenance Program	Optional	
BMS Connectivity	N/A	
Port Entry	N/A	

No alterations or additions can be made to this cabinet other than what is listed as optional