

END OF FINANCIAL YEAR EQUIPMENT *Sale*

ONE-OFF REDUCED PRICING
WHILE STOCKS LAST.



Ensure you get the best possible tax return with The Small Business Support – **\$20,000 instant asset write-off measure** announced by the Australian Government in place until June 30, 2025. **Take advantage while you still can.**

An EBOS Healthcare Publication
Primary Care Division Phone: 1800 269 534
E&OE. Offers not available to medical wholesalers and
brokers. All offers valid until: 30th June 2025 unless otherwise
specified or while stocks last.

ALL ENCLOSED PRICES ARE GST EXCLUSIVE

 INDICATES SHIPPING COSTS ARE ADDITIONAL

 **EBOS**
HEALTHCARE

ESSENTIAL TO HEALTH

- 10 % GST will apply to all items unless they are marked with GST FREE
- Freight charges apply for orders under \$250 (GST excluded)
- Prices shown are valid from April until 30th June 2025

MESI mTABLET ECG Diagnostic System

- The most powerful and simultaneous 12-lead acquisition with the University of Glasgow ECG Analysis Algorithm
- Completely wireless and very portable diagnostic device
- It captures and analyses high-quality ECG data in real time that can be shared and accessed for consultation from any device.
- Advanced diagnostic tools and analysis options available on MESI mRECORDS web platform
- Large screen offering portrait and landscape view with various lead layouts



\$5,250

EBOS Code: 11000476

MESI mTABLET ABI Diagnostic System

- PADsense™ algorithm for detection of severe Peripheral Arterial Disease
- Quick and reliable ABI measurement with pulse waveform interpretation



\$6,250

EBOS Code: 11000475

MESI mTABLET Spirometry Diagnostic System

- Multiple measurement modes and parameter calculation options
- Automatic selection of the best of all repeated maneuvers with the BestBreath™ detection
- Detailed report with the ability to switch between charts and values for a clear interpretation
- Option to compare different measurements on the same screen through MESI mRECORDS
- Pneumotachograph technology with integrated self-calibration for accurate measurements at any time

*Transducers not included with mTABLET Spirometry Diagnostic System



\$4,950

EBOS Code: 11000441





MESI mTablet ECG Module & 4-Port Charging Station

- Glasgow Interpretation Algorithm included
- ECG signal filters (MESI signal enhancement, high-pass, low-pass, mains, myogram)
- Recording speed and sensitivity settings
- Customisable printouts directly from the mTABLET
- Award-winning design and user-friendly interface
- Patient history and trends already on the tablet



\$3,850

EBOS Code: 11000500

MESI mTABLET ABI Diagnostic Module

- PADsense™ algorithm for detection of severe Peripheral Arterial Disease.
- Quick and reliable ABI measurement with pulse waveform interpretation.
- Pulse waveforms and oscilation graphs
- Advanced review and alert



\$5,250

EBOS Code: 11000538

MESI Spirometry Extension Module for mTABLET, Wireless

The MESI Spiro Extension Module allows general practitioners to perform spirometry with an existing MESI mTablet system, offering essential patient monitoring in a wireless manner. The digital Spirometer utilises Pneumotachograph technology with integrated self-calibration for accurate measurements at any time.



\$3,250

EBOS Code: 10010247



- 10 % GST will apply to all items unless they are marked with GST FREE
- Freight charges apply for orders under \$250 (GST excluded)
- Prices shown are valid from April until 30th June 2025

Examination Couch - 3 Section

Specifications:

Length	
Total	195cm
Head Section	63cm
Middle Section	63cm
Leg Section	63cm
Width	71cm
Height Adjustable	46.5cm - 92.5cm
Loading Weight	225kg
Padding thickness	5cm
Foam Density	36kg/m3



EBOS Code	DESCRIPTION	UOM	PRICE
10016248	3 SECTION ELECTRIC COUCH NAVY	Each	\$1,199
10016247	3 SECTION ELECTRIC CHARCOAL GREY	Each	\$1,199



Three Section Medical Treatment Couch Hi Lo Electric Grey

- Three Equal Sections
- Head Section with Face Hole and Plug
- Upholstered with heavy-duty PVC vinyl and supportive Foam Padding
- Rounded frame legs, adjustable feet for uneven surfaces, and 7.5cm retractable castors
- Reliable TiMotion Motor



\$1,199

EBOS Code: 10022778



FREE
Bed Roll
Holder Valued
at \$70

2 Section Hi Lo Electric Exam Couch

- Electric Hi Lo Exam Couch operated by foot switch
- Sleek design with rounded frame legs for safety
- Minimum height 50cm to maximum height 90cm in 30 seconds
- Adjustable head section through dual gas operated struts from -20° to +80°
- Head section has face hole and plug
- Adjustable feet for uneven surfaces
- 7.5cm castors for safety and mobility
- Retractable dual locking castors on all 4 corners



\$1,099

EBOS Code: 10016243



FREE
Bed Roll
Holder Valued
at \$70

Emergency Trolley

- Defibrillator shelf
- Extension table
- 304 Stainless steel IV pole
- Cardiac CPR board
- Waste bin
- External power holder & plug board
- Oxygen tank holder
- 304 Stainless steel wire basket
- Central or seals lock
- Castors with brakes
- Non slip mat
- Luxurious push handle
- Soft closing drawers
- Dimensions: 650 x 480 x 980mm
- Drawer Sizes:
3 x 110mm, 1 x 135mm, 1 x 195mm
- G.W 55kg
- Central locking system through side key



ALL ACCESSORIES INCLUDED

\$1,250

EBOS Code: 10011835



Powder Coated Trolleys

- 50 x 50 x 96.7 cm



EBOS Code	DESCRIPTION	UOM	PRICE
10011845	POWDER COATED TROLLEY 1 DRAWER	Each	\$310
10011846	POWDER COATED TROLLEY 2 DRAWERS	Each	\$350

Stainless Steel Trolleys

- 50 x 50 x 90 cm



EBOS Code	DESCRIPTION	UOM	PRICE
10011855	STAINLESS STEEL TROLLEY 1 DRAWER	Each	\$425
10011857	STAINLESS STEEL TROLLEY 2 DRAWER	Each	\$550



Welch Allyn DuraShock DS66 Trigger Aneroid Sphygmomanometer with Multi-Cuff

- Gear Free
- Shock Resistant
- Jewel Movement
- 15-Year Calibration Warranty
- Zipper Case
- Four-cuff Kit: Size-09 (Child Print) to Size-12
- FlexiPort Reusable
- 1-Tube cuffs



\$285

EBOS Code: 10009053

Welch Allyn ProBP™ 2000 Digital Blood Pressure Device

- Blood pressure readings are captured on inflation in about 20 seconds
- Meets AAMI standards for blood pressure accuracy
- Compatible with Welch Allyn FlexiPort Reusable Blood Pressure Cuffs, from child size 9 up to thigh size 13 (15 - 55cm range)
- Device ships with one adult size 11 FlexiPort Reusable Blood Pressure Cuff



\$165

EBOS Code: 10009157

Braun Thermoscan® PRO 6000 Ear Thermometer with Small Cradle

The Braun Thermoscan® PRO 6000 Ear Thermometer featuring advanced measurement technology, has the speed and simplicity you need to capture accurate and reliable temperature readings through our innovative PerfectTemp™ technology.



\$199

EBOS Code: 10008847

Welch Allyn 3.5V Portable Diagnostic Set, Coaxial LED Hard Case

The Welch Allyn 3.5V Diagnostic Set includes a coaxial LED ophthalmoscope, MacroView Basic LED Otoscope, one Lithium-Ion Rechargeable power handle and a hard case. The Diagnostic set helps keep the tools physicians need fast, organized and at hand to help examine the eyes, ears, nose and throat. The set features long-lasting power handles and LEDs, for more dependable tools to keep physicians running.



\$695

EBOS Code: 10009023

Green Series 300 General Exam Light with Mobile Stand

Welch Allyn® Green Series™ 300 Exam Light Utilising a single white LED, the GS 300 Exam Light provides a whiter, brighter light than halogen and consumes less energy. It features a perfect combination of intense light output and cool operation in a compact design, along with touchless switch to reduce risk of cross-contamination, and multiple mounting operations for floor space maximisation. It has a low cost of ownership due to elimination of bulb replacement and lower energy consumption, and easy assembly.



\$695

EBOS Code: 10009023

High Frequency Desiccator Bovie DERM 942

- 40 watt frequency desiccator
- Adjustable in 1/10th watt increments in the first 10 watts, in monopolar and bipolar
- More accurate power with Universal Power Control 100-240 VAC
- Less to do and less to buy makes this the most cost effective and user friendly unit on the market
- The Bovie DERM 942 is a ground reference unit which means no grounding pad is required



\$1,695

EBOS Code: 10016686

- 10 % GST will apply to all items unless they are marked with GST FREE
- Freight charges apply for orders under \$250 (GST excluded)
- Prices shown are valid from April until 30th June 2025

Nuline PR350 Vaccine Refrigerator 350L

- Fan Forced System
- Top Mounted Door Lock
- Interior LED
- Digital Display with Built-in Alarms
- Double Glazed Door
- Hi/Lo Readout
- Door Fan Switch
- Highly Energy Efficient
- QCPP Compliant



\$2,550

EBOS Code: 11000665

Nuline HR200G Pharmacy Refrigerator 135L

- QCPP Approved
- Fan Forced System
- Self Closing Glass Door
- Door Lock
- Interior LED Lighting
- Data Logger included
- Digital Display with Built-in Alarms
- Hi/Lo Read Out
- Remote Alarm Contacts
- Door Fan Switch
- Battery Back Up for Controller
- Port Entry
- 4 Shelves, 3 Large and 1 Small



\$1,750

EBOS Code: 11000472

Nuline Secure Pharmacy Vaccine Fridge HR400G With Glass Door 350L

The Nuline Pharmacy Vaccine Fridge is designed with an automatic defroster along with an in built alarm system for temperature faults. The Vaccine Fridge uses both LED lighting strips and double glazed glass door for energy efficiency. It is constructed to be secure with a self closing and lockable door and is suitable for pharmaceutical storage.

- Dimensions: 595 x 650 x 1850cm



\$2,850

EBOS Code: 11000738

Seca 803 Scale

Seca Clara 803 Flat Digital Scale boasts a two-component high-quality rubber coating for comfortable use as well as four stable rubber feet to safely stand on the scale. The 28-millimeter LCD display shows large figures with its functions of a tap-on automatic switch-on and an automatic switch-off. The scale also displays its figures in kg or lbs. The scale has a weight capacity of 150kg (330 lbs) as well as a graduation weight of 100g.



\$99

EBOS Code: 10011810

Seca 813 Scale

The Seca 813 Electronic Floor Scale can handle up to 200 kilograms on its extra-wide platform which makes this suitable for bariatric patients.

The tough rubber coating of the Seca 813 scale has a high capacity and the flat, low construction makes it easy to step on. It features a large LCD display for easy reading. Kg/lbs switch-over functions with 2 year warranty.




\$195

EBOS Code: 10011813

seca Digital Column Scale 250kg with Height Rod

With eye-level display is the medical-grade solution for measuring and weighing. It features a tilt-proof full metal platform and robust single piece column, that won't give way when pushing the keypad.



\$730

EBOS Code: 10011082



Delta 30 Pro Dermatoscope Bundle - Small

DELTA 30 PRO bundle with contact plate with scale, small contact plate, E4-USBC power supply with cable, charger 30 table charging station, universal smartphone adapter, case.



\$2,450

EBOS Code: 10021780

Defibtech Lifeline Semi Automatic Two Button AED Defibrillator

The Defibtech Lifeline AED is an award-winning, lightweight, and portable device designed with maximum visibility in mind. The Lifeline features a simple two-button interface, clear voice prompts, and easy-to-see progress lights.



\$1,195

EBOS Code: 11000525

Defibtech Lifeline VIEW Defibrillator

The Fully Automatic Lifeline AED is designed to analyse heart rhythms and take action, automatically delivering a shock, without the need for intervention.



\$1,195

EBOS Code: 11000527

Welch Allyn® Desk Set PanOptic MacroView Basic Ophthalmoscope

Compact, portable unit that provides convenient drop-in recharging of Welch Allyn desk/well-charger instrument power handles.



\$895

EBOS Code: 10008906

Green Series 900 Procedure Light with Ceiling Mount

- Six White LEDs
- 80,000 LUX @ 24in/61cm
- 5,500° K
- Touchless Controls
- 100-240V 50-60 Hz AC



\$1,995

EBOS Code: 10009116

Walmay Cryospray Gun 300ml

Primary care physicians can simply and quickly treat such problems as keratoses, condylomata, warts, skin tags and other skin lesions every day using Walmay Cryospray. Not only is the procedure simple, there are many benefits for both the physician and the patient in its use.



FREE
Skin Protector
with Every
Flask

\$775

EBOS Code: 10008862

Dewar Liquid Nitrogen Container, 20 Litre

The Dewar Liquid Nitrogen Container is made of a high-strength aluminium alloy with advanced powder coating technology and designed for low liquid nitrogen consumption. Each Dewar Container features a complete safety lock lid, neck plug, leather (PU) bag & ladle.



FREE
Liquid Nitrogen
Extractor Valued
at \$200

\$995

EBOS Code: 10009495

- 10 % GST will apply to all items unless they are marked with GST FREE
- Freight charges apply for orders under \$250 (GST excluded)
- Prices shown are valid from April until 30th June 2025



SPECIAL PRICING ON STETHOSCOPES

**Great savings
on 3M Littmann
Stethoscopes**



Classic II Infant Stethoscope Latex-Free Black 27 inch

- Chestpiece: Dual-sided with a 2.7 cm diaphragm
- Eartips: Soft-sealing, snap-tight, with optional large firm eartips
- Tubing: Next-generation, latex-free
- Length: 27 inches individual fit and comfort. Angled eartubes align with ear canals.



\$149⁹⁹

EBOS Code: 10021105

Cardiology IV Hunter Green Stethoscope Latex Free 69cm

- Tube Colour: Hunter Green
- Chestpiece Finish: Machined Stainless Steel
- Length: 69 cm
- Chestpiece Weight: 87g
- Diaphragm Diameter: 4.3cm



\$299⁹⁹

EBOS Code: 10020995

Classic II Paediatric Stethoscope Caribbean Blue

- Two sided chestpiece
- Colour: Caribbean Blue
- Length: 28 inches
- Diaphragm Diameter: 3.3cm with ear canals.



\$149⁹⁹

EBOS Code: 10021106

Stethoscope Classic III Stethoscope Latex Free Black/Champagne

- Colour: Black/Champagne
- Length: 69cm
- Chestpiece Weight: 82g
- Diaphragm Diameter: 4.3cm



\$179⁹⁹

EBOS Code: 10020987

ORDISI LED Examination and Treatment Room Lighting

With a small light head with integrated handle it is easy to move without blocking the vision of the operator. The flexible goose neck allows positioning in a diverse range of angles for optimal illumination of the work area.



EBOS Code	Description	UOM	PRICE
10012208	ORDISI MOBILE LED EXAMINATION AND TREATMENT LIGHT	Each	\$425
10012200	ORDISI WALL MOUNTED EXAMINATION LIGHT WITH SWIVEL BRACKET	Each	\$375

ORDISI Procedure Lamp IGLUX Series

The ORDISI Iglux procedure lamps offer higher colour rendering (CRI) and minimize heat emissions thereby improving performance in the surgical field.



EBOS Code	Description	UOM	PRICE
10012206	ORDISI PROCEDURE LAMP IGLUX SERIES MOBILE	Each	\$2,650
10012202	LIGHT ORDISI IGLUX CEILING 80CM POLE	Each	\$2,850

Hurry!

**TAKE ADVANTAGE OF ONE-OFF
REDUCED PRICING WHILE
STOCKS LAST.**