

# MICRO-TOUCH<sup>®</sup> NITRA-TEX<sup>™</sup> E.P.

Performance Series

## Nitrile, powder-free, extended cuff examination glove with strong barrier protection

This non-sterile glove is made from a softer synthetic formulation to offer a comfortable fit and feel. It also has an excellent resistance from a broad range of commonly used chemicals. High tear resistance ensures optimal protection, and a chlorinated inner coating facilitates ease of donning. With an extended cuff, this glove is suitable for any situation where additional wrist protection is required.

### KEY FEATURES AND BENEFITS:

- Additional length for extra protection.
- Excellent resistance to a broad range of commonly used chemicals
- Protects from latex Type I allergy
- Tested for use with chemotherapy drugs\*

\* See product packaging or contact Ansell Customer Service for specific chemotherapy drug permeation times and recommendations

**WARNING:** No glove completely protects against all chemicals. Users should test the suitability of this product against the specific chemicals and environment where used.

### Recommended For

- Extra protection over the wrist and arm
- Sustained, high-risk examination procedures
- Environmentally friendly facilities
- Tested for use with chemotherapy drugs
- Protection from Type I latex allergy in HCW's or patients

### FEATURED TECHNOLOGIES



## MICRO-TOUCH<sup>®</sup> Everyday, Everywhere Protection

MICRO-TOUCH<sup>®</sup> premium medical examination gloves are scientifically designed to complement every need of medical professionals. Sophisticated design features include Ansell ergonomic technologies that ensure comfort during extended procedures, wet grip innovations that prevent harmful mishandling, and micro-thin films for maximum sensitivity. Caregivers who are fully prepared and comfortably protected by versatile MICRO-TOUCH styles have the peace of mind to focus on what matters most: the health and protection of themselves and their patients.

# MICRO-TOUCH<sup>®</sup> NITRA-TEX<sup>™</sup> E.P.

Performance Series

## PRODUCT INFORMATION

Material	Nitrile	Not Made From Natural Rubber Latex	Yes
Color	Blue	Cuff Length	Extended
Powder Content	Powder-Free	External Glove Surface	Fully Textured
Freedom from Holes (Inspection level I)	1.5 AQL	Palm Thickness (mm/mil)	0.11 / 4.3
Finger Thickness (mm/mil)	0.16 / 6.3	Allergy Prevention	Latex (Type I)
Available Sizes	XS (5.5 - 6), S (6.5 - 7), M (7.5 - 8), L (8.5 - 9), XL (9.5 - 10)	Tested For Use With Chemotherapy Drugs	Yes, in accordance with ASTM D6978 and US FDA cleared
Sterile	No	Antistatic	Not Tested
Glove Length (mm/inches)	295 / 11.6	Product Segmentation	Performance Series

## ORDERING INFORMATION

Size	XS (5.5 - 6)	S (6.5 - 7)	M (7.5 - 8)	L (8.5 - 9)	XL (9.5 - 10)
Product Code	4699	4690	4691	4692	4693

## PACKAGING AND STORAGE

Packaging	100 gloves per dispenser 10 dispensers per case 1000 gloves per case
Shelf Life	3 Years
Storage Instructions	Keep out of direct sunlight; store in a cool and dry place. Keep away from sources of ozone or ignition.

## STANDARDS AND CERTIFICATIONS

ASTM D6319, Category III, EN 374:2003, EN 420:2003 + A1:2009, EN 421:2010, EN 455 1-4, EN ISO 374-1:2016, EN ISO 374-5:2016, ISO 10993-10:2002, ISO 9001, ISO 11193-1:2008, ISO 13485, Food Type: Rubber

Contact your Ansell representative for ordering or more information.

North America  
US Tel: 800 952 9916  
CA Tel: 1-844-494-7854  
Email: insidesalesus@ansell.com

Latin America & Caribbean  
Tel: +52(442) 296 20 50  
Email: cslac@ansell.com  
Brazil Tel: +55-11-3356-3100  
Brazil Email: luvas.medicas@ansell.com

Europe, Middle East & Africa  
Tel: +32 (0) 2 528 74 00  
Email: info@ansell.eu

Asia Pacific  
Tel: +603 8310 6688  
Email: apac.medical@ansell.com

Australia & New Zealand  
Tel: +61 3 9270 7270  
Email: protection@ap.ansell.com

Ansell, <sup>®</sup> and <sup>™</sup> are trademarks owned by Ansell Limited or one of its affiliates. US Patented and US and non-US Patents Pending: [www.ansell.com/patentmarking](http://www.ansell.com/patentmarking).  
© 2021 Ansell Limited. All Rights Reserved.