



## Alaris™ GP Plus Series primary set

63420EB



- Needle-Free Connectors
- Set Length 285 cm
- Priming Volume 26 ml
- Bore 3 mm
- 1 Spiked Drip Chamber(s)  
(Vented/Unvented)
- Hydrophobic Air Inlet Filter (with Cap)
- Drops per ml 20
- 2 SmartSite® Valve Y Port
- 1 Roller Clamp(s)
- 1 Slide Clamp(s)
- Pump Segment
- Alaris™ Safety Clamp
- 1 Back-Check Valve(s)
- 1 Male Luer(s)
- Rotating Male Luer(s)
- Natural rubber latex is not part of the material formulation
- DEHP is not part of the material formulation
- Quantity 50



50



20



2 SmartSite®  
Valve Y Port



1 Back-Check  
Valve(s)



1 Male Luer(s)



0086