

DISPOSABLE TOURNIQUET LATEX FREE

For improved infection control

Haines®
MEDICAL AUSTRALIA

Haines® Medical Australia offers a simple, cost effective alternative to reusable tourniquets. Many of our customers are concerned that their reusable tourniquets are:

- not being assessed properly after each use
- when soiled, still being reused posing an infection risk to the next patient
- old / soiled presenting a poor image of the facility

WHY SWITCH TO DISPOSABLE TOURNIQUETS?

A recent study examined 100 reusable tourniquets after use on patients' arms and found that:

"61 per cent of them were colonised with bacterial species not normally found in the upper arm area. Multi-resistant organisms (MRO), including the superbug methicillin-resistant *Staphylococcus aureus* (MRSA), were found in one quarter of the tourniquets."

Source: Medical Journal of Australia

USED TO

- reduce infection risk between patients
- reduce costs and time associated with cleaning of reusable tourniquets
- improve access to a fresh tourniquet (box of 100 ready to use)



Easy to use dispenser box for each 100 Tourniquets

Disposable Tourniquet

- 2.5cm x 50cm x 0.5mm thick
- Latex Free
- Soft and comfortable for patients
- Excellent fit on all limb shapes and sizes due to its stretch and length
- Made from Latex Free, Thermoplastic Elastomers (TPE)
- Flame Retardant - AS 2755.2-1985

Code: **DT100**

BOX 100 **\$25** (25c ea)
GST FREE

Code: **DT2000** (Box 2000 - 20 Boxes of 100)

BOX 2000 **\$390** (19.5c ea)
GST FREE



If you are working for a company or organisation within healthcare, aged care or a hospital, free samples are available. To request a free sample contact Haines® Medical



3,3 PRICES CURRENT FROM 1ST JULY 2018

HAINES® MEDICAL AUSTRALIA. ALWAYS THINKING. ALWAYS THERE.

Specialists in Disposable and Single Patient Use products

26 Heath Street, Lonsdale
South Australia 5160

T 08 8294 5999
F 08 8294 4337

E sales@hainesmedical.com.au
W hainesmedical.com.au

Tourniquets

