

# SEVILLE

## MC310 - FOLDING FLOORLINE BED

The SEVILLE offers the user and carer an easy to install and transportable floorline bed which is ideal for use in most nursing environments.

Features include:

- Four section mattress frame
- Integrated 200mm bed extension
- Lockout handset
- Battery backup
- Complete trendelenburg functions



Excellent under bed clearance

Easy to clean plastic / laminate surfaces

Removable plastic mattress platform covers

I.V / Self help pole sockets



Folded bed dimensions  
W: 920mm  
D: 460mm  
H: 1160mm



### Bed specifications:

**Maximum user weight**  
225kg

**Dimensions**  
L: 2320mm (2520mm extended)  
W: 920mm  
H: 1070mm

**Effective lift range**  
130 - 725mm

**Trendelenburg range**  
- 15° to + 15°

**Backrest angle**  
0 - 65°

**Leg break angle**  
0 - 30°

**Calf support angle**  
0 - 30°

**Movement**  
Linak actuation  
(including battery backup)

**Castors**  
Ø100mm total lock castors

**Finish**  
Mild steel framework  
Epoxy powder coat  
ABS covers

**Bed weight**  
120kg

**Mycomfort**

A: 100 Mulgul Road, Malaga, Western Australia 6090  
T: +61 8 9248 4444 F: +61 8 9248 4450  
Consumer Help line: 1800 998 430  
E: customerservice@kcare.com.au  
W: www.kcare.com.au  
ABN: 70 161 963 257  
ACN: 161 963 257

**KCare**  
Healthcare Equipment