



Reduces
Oedema



Reusable



Cost
Effective

Lymphlex Short Stretch Compression Bandage is constructed from 100% cotton and designed to provide a sub-bandage pressure of 30-40mm Hg when applied at full stretch, with 50-70% overlaps.

Lymphlex Short Stretch Compression Bandage is an individually replenishable component of the Lymphlex Reduce Compression Therapy Kit.

COST EFFECTIVE

Washable and reusable up to 50 times*

PROVIDES THERAPEUTIC PRESSURE

Ideal for use in phlebological conditions, lymphoedema management and soft tissue injury support

30-40MM HG SUB-PRESSURE

Provides a sub-bandage pressure of 30-40mm Hg when applied at full stretch, with 50-70% overlaps



Suitable for:

Lymphoedema management

Venous leg ulcer and oedema management

Reorder	Size	Sterility	Pcs/Unit	Units/Ctn
LOR003	8cm x 5m	Non-Sterile	1	300
LOR004	10cm x 5cm	Non-Sterile	1	240
LOR005	12cm x 5cm	Non-Sterile	1	120

*Independent laundering testing Australian Wool and Textile Association

