



Conforms to Body



100% Cotton



Medium Weight

Sterile and medium weight construction, the Medicrepe (All Cotton Crepe Bandage) is made of quality fabric offering support without obstructing muscular or joint mobility.

### STERILE

Suitable for intra-operative and post-operative bandaging

### CONFORMS TO BODY CONTOURS

Allowing flexibility without restrictions of movement

### DURABLE AND HIGH STRETCH

Provides ideal padding/protection and support



### Suitable for:

Post-operative use

Retention bandaging

Securing dressings

| Reorder | Size         | Features    | Sterility | Pcs/Unit | Units/Ctn |
|---------|--------------|-------------|-----------|----------|-----------|
| ACS001  | 2.5cm x 1.6m | Unstretched | Sterile   | 150      | 1         |
| ACS002  | 5cm x 1.6m   | Unstretched | Sterile   | 120      | 1         |
| ACS003  | 7.5cm x 1.6m | Unstretched | Sterile   | 80       | 1         |
| ACS004  | 10cm x 1.6m  | Unstretched | Sterile   | 70       | 1         |
| ACS005  | 15cm x 1.6m  | Unstretched | Sterile   | 50       | 1         |