



Waterproof
Dressing



Moisture Vapour
Permeable



Transparent
Dressing

Osmocol Hydrocolloid Dressings provide optimal wound healing conditions while allowing you to monitor the healing progress – reducing the need for unnecessary dressing changes.

TRANSPARENT DRESSING

Monitor healing progress without needing to disturb or remove dressing

PROTECTIVE FILM LAYER

Bacteria and waterproof dressing

GEL FORMING

Forms gel when in contact with exudate and protects newly formed tissue



Suitable for:

Lightly to moderately exuding wounds

Prevention of superficial skin damage and treatment of wounds

Venous leg ulcers

Skin donor and surgical sites

Pressure injuries and partial thickness burns

| Reorder | Size | Pcs/Unit | Units/Ctn |
|-------------|-------------|----------|-----------|
| THIN | | | |
| AWH001 | 5cm x 5cm | 10 | 10 |
| AWH002 | 5cm x 10cm | 10 | 10 |
| AWH003 | 5cm x 20cm | 10 | 10 |
| AWH007 | 5cm x 25cm | 10 | 10 |
| AWH004 | 10cm x 10cm | 10 | 10 |
| AWH005 | 15cm x 15cm | 5 | 10 |
| AWH006 | 20cm x 20cm | 5 | 10 |

| Reorder | Size | Pcs/Unit | Units/Ctn |
|------------------------|-------------|----------|-----------|
| STANDARD | | | |
| AWH050 | 10cm x 10cm | 10 | 10 |
| AWH051 | 15cm x 15cm | 5 | 10 |
| AWH052 | 20cm x 20cm | 5 | 10 |
| BORDERED | | | |
| AWH075 | 10cm x 10cm | 10 | 10 |
| AWH100 | 14cm x 14cm | 5 | 10 |
| SACRUM BORDERED | | | |
| AWH150 | 15cm x 18cm | 5 | 10 |