

3 Way Stop Cock



PRODUCT DATA SHEET

DESCRIPTION

The 3-Way Stop Cock allows for 360° rotation with a transparent housing for greater visibility. It requires only a minimum amount of priming fluid and can take pressures up to 4.48bar (65Psi). The arrows on the handle indicate the direction of flow.

FEATURES AND BENEFITS


- 360° rotation
- Transparent channel
- Minimum priming volume required
- Male luer lock with rotating collar
- Two female luer connectors
- Pressure up to 4.48bar (65Psi)
- Arrow on handle indicates direction of flow
- DEHP free

PRODUCT DATA

- Priming volume: 0.3mL
- Shelf life: 5 Years
- Material: Polycarbonate
- ARTG Number: 278076
- Colour: White
- Latex free
- PVC free
- Single use only
- The product is MRI and CT suitable

PRODUCT PACKAGING LABEL

RE-ORDER
IV110000



DESCRIPTION
3-Way Stop Cock
White

LOT **88888888**
yyyy/mm
yyyy/mm
00 kg

QUANTITY
500 Pieces/Ctn (10 boxes, 50pcs/box)

STERILE

m|devices
 Multigate Medical Devices Pty Ltd
 27 Llewellyn Avenue Villawood NSW 2163 Australia
 T +61(0)2 8718 2800 F +61(0)2 8718 2801
 E info@mdevices.com.au
 www.mdevices.com.au

POLYMERBOND LTD.
 100, SOUTH STREET, WILLOWOOD, NSW
 AUSTRALIA
 15000000000000000000000000000000

RE-ORDER	SIZE	CARTON QTY
IV110000		50