



Strong and Durable



Various Sizes



Colour Change Indicator

Produced with medical grade paper on one side and a PET/PP blend film on the other, Bolsaplast Sterilisation Flat Reels offer extremely low fibre shedding and are resistant to film tearing. The reels are manufactured under strict quality control guidelines to ensure consistent package integrity and indicator change after each sterilisation cycle.

COLOUR CHANGE INDICATOR

Chemical indicators for Steam/EO/Formaldehyde

NON-TOXIC INK

No damage to instruments

MEDICAL GRADE CONSTRUCTION

Offers excellent peelability and is resistant to film tearing



Suitable for:

Sterilising devices and equipment

| Reorder | Size | Pcs/Unit | Units/Ctn |
|------------|---------------|----------|-----------|
| MRP1040 | 5cm x 200m | 1 | 6 |
| MRP1041 | 7.5cm x 200m | 1 | 4 |
| MRP1042 | 10cm x 200m | 1 | 4 |
| MRP1043 << | 12.5cm x 200m | 1 | 4 |
| MRP1044 | 15cm x 200m | 1 | 4 |
| MRP1045 | 20cm x 200m | 1 | 2 |
| MRP1046 | 25cm x 200m | 1 | 2 |
| MRP1047 | 30cm x 200m | 1 | 2 |
| MRP1048 | 40cm x 200m | 1 | 1 |

<< Indent Item