



The most trusted bacterial and viral filter
on the Australian and New Zealand
markets for over 2 decades



BIRD HEALTHCARE

www.birdhealthcare.com

Key Features:

- *Manufactured from inert, latex free, eco-friendly polypropylene plastic*
- *Utilises world leading filter material for superior filtration at high flow rates*
- *All in one integrated elliptical mouthpiece provides maximum seal, comfort and prevention of leaks*
- *Colour coded for easy identification between multiple systems*

Material *GVS (formerly Airsafety)*

Filter Media *TEC 200gm Electrostatic filter media. Electrostatic hydrophobic double outer scrim*

Housing *Polypropylene plastic (100% recyclable)*

Mouthpiece Shape *Oval, Elliptical*

Mouthpiece Adapters *Paediatric/small mouthpiece adapters available*

Dead space *50ml*

Resistance *0.71cmH₂O/l/s at a flow of 12L/s*

Filter efficiency *Bacterial/Viral filtration 99.9999%*

Compatible with *Jaeger, Sensormedics, Medisoft, Ganshorn, MedGraphics, nSpire, Ndd, Vitalograph, Thorasys, CosMed, MIR.*

Packaging *100 per box (Grip nose clips optional)*

Patented Ergonomic Mouthpiece



Paediatric Mouthpiece Adapter Also Available



Grip Nose Clip Also Available



The SureGard® filter is an Australian owned, highly effective filter barrier device. The SureGard® provides both technicians and patients with optimal protection from cross contamination and infection. The R&D of the SureGard® is attributed to contributions from the medical and scientific community of Australia and New Zealand.



RJVKB1 (Pink)

Direct fit nSpire, MedGraphics (internal fit), MicroMedical, MIR SpiroBank, MIR MiniSpir, MIR Spiro Lab, MIR SpiroDoc Fukuda and Peak flow meters.
External Diameter- 28.5mm
Internal Diameter- 25mm



RJVKB2 (Purple)

Direct fit Viasys, SensorMedics, Jeager, QRS SpiroCard.
External Diameter- Slight taper from 33mm to 32mm
Internal Diameter- 29mm



RJVKB3 (Blue)

Direct fit MedGraphics (internal fit), MicroMedical, MIR Spiro Lab, MIR SpiroDoc, MIR SpiroBank, MIR MiniSpir, Peak flow meters, Viasys, SensorMedics, Jeager, QRS SpiroCard.
External Diameter (right)- Slight Taper from 33mm to 32mm
Internal Diameter (right)- 29mm
External Diameter (left)- 28.5mm
Internal Diameter (left)- 25mm



RJVKB4 (Green)

Direct fit Vitalograph Spirometers, MicroMedical, Spider USb, MicroLoop, SpiderXpert, MIR systems & CareFusion
External Diameter- 32.5mm
Internal Diameter- 29mm



RJVKB5 (Red)

Direct fit for all MediSoft systems
External diameter- 31mm
Internal diameter- 29mm



RJVKB6 (Clear)

Direct fit Ventilation systems, nebulisers, CPAP & BiPAP equipment.
External Diameter (right)- 25.5mm
Internal Diameter (right)- 22mm
External Diameter (left)- 22mm
Internal Diameter (left)- 19.5mm



RJVKB7 (Pearl)

Direct fit ndd Easy One Spirometer (spirettes)
External Diameter (right)- 20mm
Internal Diameter (right)- 18.5mm



RJVKB8 (Orange)

Direct fit for all MedGraphics systems (external fit)
External diameter- 36mm
Internal diameter- 34mm



RJVKB9 (White)

Direct fit for all CosMed Pony, CosMed MicroQuark
External diameter- 30mm
Internal diameter- 28mm



RJVKB10 (Yellow)

Direct fit nSpire KoKo
External diameter- 47.5mm
Internal diameter- 45mm




RJVKB11 (Orange/White)

Direct fit TremeFlo FOT
External diameter- 36mm
Internal diameter- 34mm



RJVKB12 (White/Orange)

Direct fit MedGraphics (Clip On)
External diameter- 36mm
Internal diameter- 34mm



BIRD HEALTHCARE

www.birdhealthcare.com

sales@birdhealthcare.com

Australia
Melbourne

New Zealand
Auckland

France
Paris

England
London

Hong Kong
Hong Kong